# PerQseal<sup>®</sup>

## CLOSURE DEVICE

En: INSTRUCTIONS FOR USE .....	3
De: GEBRAUCHSANWEISUNG.....	12
Fr: MODE D'EMPLOI .....	22
It: ISTRUZIONI PER L'USO .....	32
Nl: GEBRUIKSAANWIJZING .....	41
Da: BRUGSANVISNING .....	50
No: BRUKSANVISNING .....	59
Sv: BRUKSANVISNING .....	68
Fi: KÄYTTÖOHJE .....	77
Es: INSTRUCCIONES DE USO .....	86
Pt: INSTRUÇÕES DE UTILIZAÇÃO .....	96

**This page is intentionally left blank**

## INSTRUCTIONS FOR USE

To ensure proper deployment and use of this device and to prevent injury to patients, read all information contained in these instructions

For assistance, contact:  
Vivasure Medical customer Service:

Telephone Number (9:00 – 17:00 GMT Mon to Fri)	+353 91 395 440
Fax Number	+353 91 395 399

Prior to use, the operator must review the Instructions for Use and be familiar with the deployment techniques associated with the use of this device.

### Warranty

VIVASURE warrants that reasonable care has been used in the design and manufacture of this device. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilisation of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond VIVASURE's control directly affect the device and the results obtained from its use. VIVASURE shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. VIVASURE neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. **VIVASURE assumes no liability with respect to devices reused, reprocessed or resterilised and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such devices.**

### TABLE OF CONTENTS

ENGLISH/ENGLISH	
WARRANTY .....	3
DEVICE DESCRIPTION .....	4
HOW SUPPLIED .....	4
INDICATIONS FOR USE .....	4
CONTRAINDICATIONS .....	4
WARNINGS .....	4
PRECAUTIONS .....	5
SPECIAL PATIENT POPULATIONS .....	5
PROCEDURE .....	5
POST PROCEDURAL MANAGEMENT .....	10
TROUBLE SHOOTING GUIDE .....	10
RECOMMENDED RECOVERY OPTIONS .....	11
RISKS .....	11
GRAPHICAL SYMBOLS FOR MEDICAL DEVICE LABELLING .....	106

### Device Description

The PerQseal® closure device is a Vascular Closure Device (VCD), designed to percutaneously deliver a bioabsorbable vascular implant to seal femoral artery puncture sites following endovascular procedures.

The implant has both intraarterial and extraarterial components. The implant is inserted into the artery through an Introducer-sheath and over a guidewire. The implant is designed to be positioned relative to the arteriotomy and effect a seal. The implant is expected to be fully absorbed within 180 days.

The PerQseal® closure device has a delivery system for the percutaneous delivery of the implant to the arteriotomy. The delivery system has a handle at its proximal end and a flexible shaft, which attaches to the implant at the distal end.

Refer to Figure 1 showing the delivery system and implant, indicating the key features of the PerQseal® closure device. The delivery system is designed to ensure single use of the device.

Figure 1. PerQseal® closure device

1. Handle	6. Release cam (Back cam)
2. Shaft	7. Loading cannula
3. Implant	8. Funnel lever
4. Sheath cam	9. Removable tray lid
5. Cam lock	

The PerQseal® closure device will be used in conjunction with the following:

- PerQseal® Introducer (Introducer-sheath and Dilator)
- A 0.014" (0.356 mm) guidewire of minimum length 150 cm – not provided
- 20 ml syringe – not provided

*Note: The term 'Device' is defined as the entire VCD system including delivery system, and implant. The term 'Product' is defined as the entire VCD system including delivery system, implant, Instructions for Use (IFU), labelling and packaging.*

### How Supplied

The packaged device is supplied in a procedure pack, in one of two models, either "S" or "L". Each procedure pack contains the PerQseal® closure device in a foil pouch, the PerQseal® Introducer (S or L) in a Tyvek®<sup>1</sup> pouch, two separate IFUs for each respective device (PerQseal® and PerQseal® Introducer) and Quick reference guide. The PerQseal® foil pouch contains a Tyvek® pouch which in turn contains the PerQseal® closure device, housed in a thermoformed tray. Note the foil pouch is a non-sterile moisture barrier. The PerQseal® Introducer pouch contains the Introducer-sheath and Dilator mounted on a support card. **The Tyvek® pouches provide the sterile barrier for their respective contents.** Contents are non-pyrogenic and supplied STERILE. All packaging not contained within the Tyvek® pouches is non-sterile. The PerQseal® closure device is sterilised using an ethylene oxide (EO) process. Please refer to the end of this booklet for graphical symbols for medical device labelling.

Store in a cool, dry, dark place.

Vivasure Medical Limited  
Parkmore Business Park West  
Galway  
Ireland.

### Indications for use

The PerQseal® closure device shall be indicated for the percutaneous sealing of a common femoral arteriotomy in patients following interventional therapeutic endovascular procedures.

The PerQseal® closure device is indicated for use with 0.014" (0.356 mm) guidewires of minimum length 150 cm only.

The PerQseal® closure device is indicated for use with either the 'S' or 'L' PerQseal® Introducer.

### Contraindications

There are no contraindications to the use of this device. Attention is drawn to the warnings, precautions, and special patient populations.

### Warnings

- Contents supplied STERILE. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilise. Reuse, reprocessing or resterilisation may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
- Do not use the Product if the packaging has been previously opened or damaged, or if any of the components appear to be damaged or defective. This may result in patient injury, illness or death.
- Discard product if 'Use By date' has passed or 'Use By date' is illegible.

<sup>1</sup> Tyvek® is a trademark of DuPont.

- Do not use the PerQseal® closure device if the puncture site has been made above the most inferior border of the inferior epigastric artery (IEA) and/or above the inguinal ligament based upon bony landmarks, as this may result in retroperitoneal bleeding or haematoma. Perform a femoral angiogram to verify location of the puncture site.
- Do not use the PerQseal® closure device if the puncture site has been made in the profunda femoris artery (PFA), superficial femoral artery (SFA) or located less than 10 mm proximal to the bifurcation of the profunda femoris and superficial femoral artery, since such puncture sites may result in incorrect implant position resulting in a pseudoaneurysm, intimal dissection, or vessel occlusion. Perform an oblique femoral angiogram to identify the bifurcation and verify location of the puncture site is within the common femoral artery.

#### Precautions

- This product is intended for use by physicians experienced in therapeutic endovascular techniques.
- Do not advance, torque or withdraw the PerQseal® closure device against excessive resistance. Excessive force used to advance, torque or withdraw the *PerQseal® closure device* should be avoided, as this may lead to significant vessel damage and/or breakage of the device.
- Remove the PerQseal® closure device from the patient if it is difficult to reposition to the arteriotomy.
- Remove the PerQseal® closure device from the patient if it is difficult to deploy.
- Do not use the PerQseal® closure device where bacterial contamination of the procedural sheath or surrounding tissues may have occurred as this may result in an infection.
- Do not use the PerQseal® closure device if the puncture is through the posterior wall or if there are multiple punctures, since such punctures may result in a retroperitoneal bleeding. Do not use the *PerQseal® closure device* if there is a vascular graft or stent within 20 mm of the access site.
- Do not modify the PerQseal® closure device.
- As with all percutaneous based procedures, infection is a possibility. Observe sterile techniques at all times when using the PerQseal® closure device. Employ appropriate groin management post procedure and post hospital discharge to prevent infection.
- If a patient has had a procedural sheath left in place for longer than 8 hours, consider the use of prophylactic antibiotics before insertion of the PerQseal® closure device.
- Do not use the PerQseal® closure device subsequent to or in conjunction (in the ipsilateral limb) with any other closure device or techniques, e.g. suturing.
- Following implantation of the PerQseal® closure device, avoid vascular access at the same location within 90 days.
- Consideration should be given to the access site: tortuosity; disease state (including calcification, plaque and thrombus); and arterial diameter relative to the largest sheath size outer-diameter to ensure a successful closure and to minimise vascular complications.
- Excessive compression may damage the implant and inadvertently cause bleeding.

#### Special Patient Populations

The safety of the PerQseal® closure device has not been established in the following patient populations:

- Patients suffering with definitive or potential coagulopathy or platelet count less than 100,000/µl.
- A measured activated clotting time (ACT) of greater than 350 seconds immediately prior to sheath removal.
- Patients that are expected to be continuously treated with anticoagulation therapy post-procedure such that their ACT reading is expected to be elevated above 350 seconds for more than 24 hours after the procedure.
- Patients with evidence of arterial diameter stenosis greater than 20 % within 20 mm of the arteriotomy.
- Patients with circumferential calcification within 20 mm of the arteriotomy.
- Patients with a concentration of fibrinogen less than 100 mg/dl.
- Patients with a known allergy to any of the materials used in the PerQseal® closure device.
- Patients judged unsuitable for surgical repair of the access site.
- Patients with a common femoral artery lumen diameter of less than 7 mm.
- Patients that have a lower extremity amputation from the ipsilateral or contralateral limb.
- Patients that have undergone a percutaneous procedure greater than 8 F in the ipsilateral leg, within the previous 90 days.
- Patients with Haematocrit less than 30%.
- Patients that have undergone a percutaneous procedure using an absorbable intravascular closure device for haemostasis, in the ipsilateral leg, within the previous 90 days.
- Patients with an acute haematoma of any size, arteriovenous fistula or pseudoaneurysm at the access site.
- Patients suffering significant blood loss/transfusion during interventional procedure or within 20 days prior to procedure requiring transfusion of greater than 4 units of blood.
- Patients with angiographic evidence of arterial laceration, dissection or stenosis within the external iliac or femoral artery before the use of the *PerQseal® closure device*.
- Patients suffering severe claudication, stenosis of the iliac artery > 50% or previous bypass surgery/stent placement in the region of the vascular access.
- Patients who are pregnant or lactating.
- Patients under the age of 18 years.

#### PROCEDURE

The PerQseal® closure device is to be used only by a licensed physician possessing adequate instruction and training by Vivasure Medical in the use of the PerQseal® closure device. Exercise care during PerQseal® closure device handling to reduce the possibility of accidental device breakage.

It is recommended to perform angiography in both the Anterior Posterior projection and Oblique projection (Right Anterior Oblique – 30 degrees from the right limb, Left Anterior Oblique – 30 degrees from the left limb) prior to ipsilateral access to identify the bifurcation and ensure that the location of the puncture will be in the common femoral artery.

It is recommended that ultrasound or angiographic guided puncture / micropuncture is used for access via the ipsilateral limb, to ensure a centre stick, away from anterior calcification, at least 10 mm proximal to the SFA/PFA bifurcation and access is within the common femoral artery.

Use clinical judgement to select the appropriate PerQseal® Introducer, “S” or “L”, based on vessel and arteriotomy size.

The *PerQseal® closure device* should be used within four hours of opening the foil pouch.

If the *PerQseal® closure device* cannot be positioned juxtaposed to the arteriotomy, due to improper orientation or patient vascular anatomy, the implant and delivery system should be withdrawn from the patient. Haemostasis can then be achieved by employing good clinical judgement.

#### Examination of Products

1. Carefully inspect the packaging to ensure there is no damage.
2. Check all Labelling to ensure the product is within its use by date.

#### Opening the packaging

1. Remove the foil pouch from the procedure pack.
2. Inspect the foil pouch to ensure that it has not been breached.
3. Tear open the top of the foil pouch from the side at the notched position.
4. Open the end of the foil pouch and remove the Tyvek® pouch.
5. Inspect the Tyvek® pouch to ensure the sterile barrier has not been breached.


6. While observing aseptic techniques, open the Tyvek® pouch by peeling apart the end with the chevron seal (the end with the 'v-shape' seal).
7. The device must be used within 4 hours of opening the foil pouch.
8. Present the exposed tray for aseptic transfer into the sterile field.
9. Have a person in the sterile field remove the tray containing the PerQseal® closure device from the Tyvek® pouch.

**Preparation**

1. Insert and position selected size PerQseal® Introducer-sheath as per PerQseal® Introducer IFU.
2. Using standard clinical practice insert a 0.014" (0.356 mm) guidewire with at least 50 cm of available guidewire protruding from the Introducer-sheath.  
*Note: If a vascular graft is present in the expected 0.014" guidewire path consider delivering the 0.014" guidewire via a pigtail catheter once the dilator has been removed from the Introducer-sheath.*
3. Ensure that the Introducer-sheath is positioned at least 4 graduations into the vessel as per PerQseal® Introducer IFU. The Introducer-sheath should be at a low angle, of less than 40° relative to the longitudinal axis of the femoral artery.(ideally <30°)
4. Perform an ACT reading to ensure that ACT is less than 350 seconds.
5. It is recommended to perform an ipsilateral site angiogram to confirm the arteriotomy puncture is within the common femoral artery and to assess presence of:
  - Dissection
  - Stenosis (greater than 20% within 20 mm of the arteriotomy)
  - Perforation other than the access arteriotomy*Note: If any of the above is confirmed by the angiogram, do not use the PerQseal® closure device.*


**Tray Loading**

- A. Whilst the tray is lying flat on the sterile preparation area, fill the recess at the end of the tray (indicated on the tray with "A") with sterile non-heparinised saline solution (approximately 20 cc) to submerge the implant, Figure 2.


**Figure 2. Fill tray recess with non-heparinised saline.**

- B. Pull off the small Removable tray lid, indicated with a "B" on the tray, by lifting the tab as shown in Figure 3 until completely removed.  
*Note: do not remove any other lids.*

**Figure 3. Remove Removable Tray Lid.**

- C. Holding the tray steady, grip the rear of the handle (release cam) and pull until a firm stop is reached (indicated on the tray with "C"), Figure 4. The implant is now folded in the loading cannula tube.

**Figure 4. Pull device handle back to firm stop.**

D. Flip the funnel lever to a vertical position to release the delivery device from the loading funnel (indicated on the tray with “D”), Figure 5.


Figure 5. Flip the funnel lever to release device from tray.

- Continue to pull back and withdraw the PerQseal® closure device from the tray.
- Visually check the integrity of the implant (within the loading cannula); there should be no creases, folds or tears, see Figure 6.

Figure 6. Inspection of Implant (within loading cannula).

- The PerQseal® closure device is now ready for use.

*Note: Use the PerQseal® closure device within 30 minutes of completing this step.*

**Procedural Steps:**

**Insert the Device**

**Load Guidewire:**

- Use non-heparinised saline to lubricate the guidewire during insertion. Insert the guidewire into the implant within the loading cannula, Figure 7 continue to insert the guidewire until it emerges from the rear of the PerQseal® closure device handle. Hold the PerQseal® Introducer-sheath at a low angle relative to the longitudinal axis of the artery (e.g. 20°) and keep the Introducer-sheath as straight as practical, to avoid kinking.


Figure 7. Loading guidewire into device cannula.

**First Lock:**

- Advance the Loading cannula over the guidewire towards the Introducer-sheath hub, Figure 8. If difficulty is experienced tracking over the guidewire, ensure the guidewire is not kinked and that it is free of coagulated blood and thrombus.  
*Note 1: Hold the guidewire at the rear of the handle during insertion to ensure the PerQseal® closure device travels over the wire and do not inadvertently insert additional wire into the patient during PerQseal® closure device insertion/advancement.*  
*Note 2: Only grip the back portion of the Loading cannula during insertion as indicated in Figure 7.*


Figure 8. Advance the PerQseal® closure device to the Introducer-sheath.

Maintaining the position of the Introducer-sheath insert the loading cannula fully into the Introducer-sheath hub. The loading cannula cap should be fully seated and locked into the Introducer-sheath hub, (with a click). The green marking being aligned confirms this step is successfully completed, Figure 9.


Figure 9. Seating Loading cannula into Introducer-sheath.

**Second Lock:**

- c. Holding the Introducer-sheath hub in one hand and the device handle in the other hand, carefully advance the device handle. Ensure the indicators on the Introducer-sheath and handle are in alignment and the loading cannula/sheath hub assembly locks (with a click) into the device handle. Note the green indicator in the window. This confirms that the PerQseal® closure device handle has been successfully locked. Figure 10.  
*Note: During insertion, maintain the guidewire stationary at the rear of the handle. Ensure the position of the sheath in the artery remains unchanged. At this stage the handle, loading cannula and sheath are all locked together as one.*


Figure 10. Insertion of PerQseal® closure device into Introducer-sheath.

Check the orientation/alignment of the PerQseal® closure device by ensuring the handle and Introducer-sheath markings are facing vertically upwards away from the patient (see Figure 11). If required, rotate the device (using the handle) to achieve the correct alignment.


Figure 11. Correct orientation of PerQseal® closure device prior to deployment.

Grip the distal end of the PerQseal® closure device handle (closest to patient). Stabilise the device in position by supporting this hand on the patient, Figure 12.

**STEP 1 Unsheathing the implant**

- a. Place hand on the sheath cam, marked '1' and rotate this grip section in the indicated direction 180° to align the dots. Ensure the hand gripping the distal end of the handle remains fixed and the handle remains stationary during this rotation, Figure 12. The implant is now exposed within the artery.

**Note 1:** the Introducer-sheath will withdraw from the patient into the handle of the PerQseal® closure device during this step to expose the implant in the artery. It is important that the handle does not advance toward the patient during rotation of the cam.

**Note 2:** If high resistance to rotation is felt prior to reaching the 180° rotation relieve the rotation pressure and try again. Do not rotate the cam backward.

**Note 3:** Do not rotate beyond 180°, this is indicated by the alignment of green dots.


Figure 12. Step 1 of PerQseal® closure device Actuation (Unsheathing the implant).

## STEP 2 Tamponade

- a. While holding the device handle at a low angle ( $\leq 30$  degrees), apply light extracorporeal pressure over the skin insertion site with the other hand. Then slowly withdraw the device from the artery, whilst watching the graduations on the Introducer-sheath, when approximately two graduations have been withdrawn, the Introducer-sheath tip will have exited the arteriotomy. Expect a reduction in tension needed to withdraw the PerQseal® closure device, after the Introducer-sheath tip exits the arteriotomy.


*Note: The digital pressure should ONLY be applied by the device operator to ensure tactile feedback occurs.*

- b. Continue to withdraw the device a further 1.5 graduations (= 3.5 graduations total). At this point, gentle tactile resistance should be felt as a result of the implant anchoring itself against the inside lumen of the artery wall. This is the tamponade position, and the implant is now controlling the bleeding from the arteriotomy. Figure 13.

*Note 1: It is important to confirm the location by gentle resistance of the implant against the arterial wall to avoid the implant being deployed within the tissue tract.*

*Note 2: Avoid excessive tension on the PerQseal® closure device. Excessive tension on the PerQseal® closure device will reduce the implant's effectiveness to achieve haemostasis or could result in the implant being pulled out of the arteriotomy.*

*Note 3: The PerQseal® closure device delivery system location at the arteriotomy may be confirmed by fluoroscopy/angiography or Doppler Ultrasonography.*


**Figure 13. Tamponade position, implant positioned juxtaposition to arteriotomy.**

- c. Slight adjustment of PerQseal® closure device position may be necessary to achieve tamponade. Ensure that an acceptable tamponade is achieved before proceeding with deployment.

**NOTE: Once an acceptable tamponade is achieved, avoid implant/sheath movement until the device is fully deployed.**

## STEP 3 Remove guidewire

- a. Locate indicator marked '3' and remove the guidewire completely from the patient and PerQseal® closure device. Guidewire removal may be done by an assistant. Figure 14.


**Figure 14. Removal of guidewire.**

*Note: The position of the PerQseal® closure device handle should remain fixed at all times. (Such that the implant is positioned at the arteriotomy with minimal tension.)*

## STEP 4. Release cam lock

- a. Locate the cam lock (marked '4') and fully depress this green button. Figure 15. Ensure it is fully depressed before proceeding to step 5.

**Figure 15. Release cam lock.**

## STEP 5 Implant Deployment and Release

- a. Place the free hand on the Release cam, marked '5' and partially rotate this grip section 90° in the direction indicated as far as the first alignment marks. Figure 16. The implant is now in the fully deployed configuration, but the implant is still attached to the device handle.  
*Note: Ensure that tamponade is still acceptable before proceeding with full deployment. Slight adjustment of PerQseal® closure device position may be necessary to achieve tamponade.*

- b. Fully rotate this grip section (marked '5') a further 90° in the direction indicated as far as the second alignment marks. This releases the implant from the delivery system.

*Note: Ensure the hand gripping the distal end of the handle remains fixed and the handle remains stationary during both stages of Step 5..*


Figure 16. Release implant.

- c. Gently withdraw the delivery system from the tissue tract. Figure 17.


Figure 17. Removal of deployed device.

- d. If suitable access is available an angiogram or duplex ultrasound may be performed immediately after *PerQseal*<sup>®</sup> closure device deployment to confirm arteriotomy sealing.

*Note: If haemostasis has not occurred, gentle manual compression may be applied proximal to the arteriotomy sufficient to control bleeding. Consult troubleshooting section.*

- e. Use standard clinical procedures to treat any skin incisions.

## Post procedural Management

Before considering patient discharge, assess for the following:

- Bleeding from the access site
- Anticoagulation thrombolytic or antiplatelet therapy
- Haematoma at the access site
- Hypotension
- Any condition requiring clinical observation
- Pain, discomfort or transitory local irritation and inflammation at the puncture site.
- Local and distal pulse deficits, ischaemia.

## Trouble Shooting Guide

### 1. **Cannot insert the loading cannula into the Introducer-sheath hub:**

Ensure that the supplied *PerQseal*<sup>®</sup> Introducer-sheath is being used and the loading cannula is not damaged.

### 2. **Difficult to advance the *PerQseal*<sup>®</sup> closure device into the Introducer-sheath:**

Ensure that the loading cannula is fully engaged with the Introducer-sheath hub, if not discard the *PerQseal*<sup>®</sup> closure device and replace with new *PerQseal*<sup>®</sup> closure device.

### 3. ***PerQseal*<sup>®</sup> closure device implant withdrawn from artery during positioning:**

Control bleeding by application of manual compression or cross-over balloon if *in situ*. Withdraw the *PerQseal*<sup>®</sup> closure device and Introducer-sheath from the patient, as per Section 1 below (Recommended Recovery Options), ensuring to leave the guidewire *in situ*. Discard the *PerQseal*<sup>®</sup> closure device and Introducer-sheath and replace with new *PerQseal*<sup>®</sup> closure device.

### 4. **Cannot achieve acceptable tamponade during implant positioning prior to deployment:**

- Avoid excessive tension on the *PerQseal*<sup>®</sup> closure device as excessive tension will reduce the implant's effectiveness to achieve haemostasis and control bleeding.
- Change the angle and/or orientation of the *PerQseal*<sup>®</sup> closure device (relative to the longitudinal axis of the artery) by using very small adjustments (approximately 5°).
- If acceptable tamponade cannot be achieved do not deploy the *PerQseal*<sup>®</sup> closure device. Follow the recover suggestion outlined in Section 1 below (Recommended Recovery Options).

### 5. **Difficult to withdraw the guidewire from the *PerQseal*<sup>®</sup> closure device:**

- Do not force.
- Ensure the guidewire is not snagged endovascularly.
- Lower the angle of the *PerQseal*<sup>®</sup> closure device (relative to the longitudinal axis of the artery).
- If it is not possible to withdraw the guidewire do not deploy the *PerQseal*<sup>®</sup> closure device and follow the steps outlined in Section 1 Recommended Recovery Options below.

### 6. **Haemostasis not achieved post deployment of the *PerQseal*<sup>®</sup> closure device:**

- Light manual compression may be applied proximal to the implant to control bleeding until haemostasis is achieved.
- Note: It is not recommended to apply manual pressure directly onto the location where implant is situated. Excessive compression may damage the implant and adversely cause bleeding.
- Ensure ACT level is below 350 seconds, use clinical judgement and reduce the effects of heparin with agents such as protamine sulphate.
- Coagulating extraarterial agents may be used.
- Temporary use of a cross-over balloon inflated directly beneath the implant may be used.
- Keep patient under observation until stable.

**7. Difficult to withdraw the device post deployment:**

The device tip may be caught on subcuticular soft tissue. Rotate the device in either direction and try to withdraw.

**8. Guidewire in situ post deployment:**

If the guidewire is still *in situ* post deployment, follow the recovery steps outlined in Section 4 Recommended Recovery Options.

## Recommended Recovery Options

1. If it becomes necessary to remove the implant from the artery whilst still attached to the delivery system then follow these recommended recovery steps:
  - i) Control bleeding by manual compression or via cross-over balloon inflation.
  - ii) Ensure Introducer-sheath is withdrawn from the artery (bring device to **Step 2 Achieve Tamponade** stage).
  - iii) Advance the *PerQseal*<sup>®</sup> closure device into the artery approximately 1 cm.
  - iv) Rotate the handle 180° so that the handle is now facing towards the patient in a downward direction.
  - v) Withdraw the handle to remove the *PerQseal*<sup>®</sup> closure device implant from the artery whilst maintaining the *PerQseal*<sup>®</sup> closure device at an approximate angle of 20° relative to the patient. Ensure the guidewire remains *in situ* during withdrawal.
  - vi) Re-cannulate the artery over the 0.014" guidewire with 4-9 F dilator.
  - vii) Exchange the 0.014" guidewire for a 0.035" guidewire.
  - viii) Open a new *PerQseal*<sup>®</sup> closure device procedure pack.
  - ix) Remove the 4-9 F dilator and reinsert a new Introducer-sheath as per *PerQseal*<sup>®</sup> Introducer IFU.
  - x) Exchange the 0.035" guidewire for a suitable 0.014" guidewire, ensuring approximately 50 cm length of guidewire is available from the rear of the Introducer-sheath.
  - xi) Follow the *PerQseal*<sup>®</sup> closure device deployment steps outlined above.
2. If the implant becomes detached from the delivery system within the artery whilst the Introducer-sheath and guidewire are *in situ* use the following recovery steps to retrieve the implant:
  - i) Withdraw the delivery system and the Introducer-sheath and ensure the guidewire remains *in situ*.
  - ii) Re-cannulate the artery over the 0.014" guidewire with a 4-9 F dilator.
  - iii) Insert a 0.035" guidewire past the implant.
  - iv) Open a new *PerQseal*<sup>®</sup> device Procedure Pack and insert Introducer as per *PerQseal*<sup>®</sup> Introducer IFU.
  - v) Remove the dilator.
  - vi) Insert a snare beyond the implant to snare the original 0.014" guidewire.
  - vii) Withdraw the snare, original 0.014" guidewire and trapped implant into the Introducer-sheath.
  - viii) Ensure the implant is within the Introducer-sheath and remove the snare, original guidewire, trapped implant and Introducer-sheath from the patient, whilst maintaining the position of the second 0.035" guidewire *in situ*.
  - ix) Insert *PerQseal*<sup>®</sup> Introducer and deploy the new *PerQseal*<sup>®</sup> closure device as per instructions above.
3. For recovery of the bioabsorbable implant from within the artery post Introducer-sheath and guidewire removal it is recommended to locate the implant via angiography. Based on the implant's location and the level of perfusion downstream of the implant, use clinical judgement to decide an appropriate strategy for recovery.
4. If the device has been fully deployed whilst the 0.014" guidewire is still *in situ* use the following recovery steps to retrieve the guidewire.
  - i) Remove the sheath and delivery system from the guidewire.
  - ii) Very slowly increase tension on the guidewire, at the same angle as deployment or approx. 30° to the artery, until movement of the guidewire occurs. Simultaneously apply light extracorporeal pressure slightly distal to the position of the implant.
  - iii) Continue to withdraw the guidewire whilst maintaining light extracorporeal pressure until fully removed.
  - iv) Inspect the guidewire to ensure the entire guidewire has been removed from the patient.
  - v) Expect some bleeding from the guidewire port in the implant. This can be managed by applying sufficient manual compression to control bleeding. Note: excessive compressing will damage the implant and adversely cause bleeding.
  - vi) Keep the patient under observation until stable.

### Risks

The following risks are the anticipated complications associated with use of the *PerQseal*<sup>®</sup> closure device:

- Bruising, oozing at the puncture site
- Transient extravasation
- Late bleed
- Retroperitoneal bleeding
- Haematoma
- Pain, discomfort or transitory local irritation and inflammation at the puncture site
- A mild local inflammatory reaction to the implant material (related to the ongoing biodegradation process) within the first 90 days. The breakdown products of the implant material are glyoxylate or glycine which is extracted via urine or as carbon dioxide respectively
- Access-site related nerve injury, vascular spasm
- Local and/or distal pulse deficits, ischemia
- Stenosis at the closure site
- Infection
- Arterial or deep vein thrombosis
- Pseudoaneurysm
- Arteriovenous fistula
- Vasovagal response

Other events that could possibly occur include:

- Occlusive intraluminal thrombus and/or emboli formation at the implant site
- Sinus formation
- Compartment syndrome
- Iliac or femoral artery dissection
- Partial or complete occlusion of the artery or peripheral arteries
- Embolisation of the intraarterial implant components
- Damage of the arterial wall (perforation or rupture)
- Allergic or adverse foreign body reaction
- Corrective intervention due to any of the above complications such as surgery and/or transfusion

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© 2018 VIVASURE or its affiliates. All rights reserved.