PerQseal[®] +

A NEW CLASS IN PERCUTANEOUS CLOSURE

Notice: Implantation of the PerQseal[®] + Closure Device should be performed only by physicians who have successfully completed Vivasure approved training

En: INSTRUCTIONS FOR USE	2
De: GEBRAUCHSANWEISUNG.....	12
Fr: MODE D'EMPLOI	22
It: ISTRUZIONI PER L'USO	32
Nl: GEBRUIKSAANWIJZING	42
Da: BRUGSANVISNING	52
No: BRUKSANSVISNING	62
Sv: BRUKSANVISNING	72
Fi: KÄYTTÖOHJE	82
Es: INSTRUCCIONES DE USO	92
Pt: INSTRUÇÕES DE UTILIZAÇÃO	102

INSTRUCTIONS FOR USE

English/English

To ensure proper deployment and use of this device and to prevent injury to patients, read all information contained in these instructions

For assistance, contact:
 Vivasure Medical customer Service:

Telephone Number (9:00 – 17:00 GMT Mon to Fri)	+353 91 395 440
Basic UDI-DI	539153318DP2FA2VCDUC

Warranty

VIVASURE warrants that reasonable care has been used in the design and manufacture of this device. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilisation of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond VIVASURE’s control directly affect the device and the results obtained from its use. VIVASURE shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. VIVASURE neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. **VIVASURE assumes no liability with respect to devices reused, reprocessed or resterilised and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such devices.**

TABLE OF CONTENTS

ENGLISH/ENGLISH.....2

Warranty.....2

Device Description.....3

How Supplied3

Intended Purpose3

Indications for use3

Clinical Benefits.....4

Contraindications.....4

Warnings.....4

Precautions4

Special Patient Populations5

Procedure.....5

Post procedural Management.....10

Recommended Recovery Options11

Adverse Events & Risks11

Device Description

The PerQseal+ closure device is a Vascular Closure Device (VCD), designed to percutaneously deliver a bioabsorbable vascular implant (with the presence of an abluminal surface matrix) to seal femoral artery puncture sites following endovascular procedures.

The implant has both intraarterial and extraarterial components. The implant is inserted into the artery through an Introducer-sheath and over a guidewire. The implant is designed to be positioned relative to the arteriotomy and effect a seal. The implant has passed design validation testing in pre-clinical and clinical studies in the presence of procedural fluids used in the handling and preparation of the implant without evidence of interaction with the absorbable material.

The PerQseal+ closure device has a delivery system for the percutaneous delivery of the implant to the arteriotomy. The delivery system has a handle at its proximal end and a flexible shaft, which attaches to the implant at the distal end.

Refer to Figure 1 showing the delivery system and implant, indicating the key features of the PerQseal+ closure device. The delivery system is designed to ensure single use of the device.

Implant Composition: <116 mg Polydioxanone. Polydioxanone is a synthetic, absorbable, low reactivity polymer, degradation of which occurs along the polymer back bone by hydrolysis of the ester bonds. This occurs randomly along the polymer chain cutting the bonds and converting the long polymer chains into (shorter fragments) of natural monomeric acids found in the body. The polymer is thus degraded to non-toxic, low molecular weight residues capable of being eliminated from the body by normal metabolic pathways. The implant begins to lose mechanical functionality after ~ 40-days with total absorption of the implant at ~ 180-days.

Figure 1. PerQseal+ closure device

1. Handle	6. Release cam (Back cam)
2. Shaft	7. Loading cannula
3. Implant	8. Funnel lever
4. Sheath cam	9. Removable tray lid
5. Cam lock	

The PerQseal+ closure device will be used in conjunction with the following:

- PerQseal Introducer (L Introducer-sheath and Dilator)
- A 0.014" (0.356 mm) guidewire of minimum length 150 cm – not provided
- 20 ml syringe – not provided

Note: The term 'Device' is defined as the entire VCD system including delivery system, and implant. The term 'Product' is defined as the entire VCD system including delivery system, implant, Instructions for Use (IFU), labelling and packaging.

How Supplied

The packaged device is supplied in a procedure pack. Each procedure pack contains the PerQseal+ Closure Device in a foil pouch, the PerQseal Introducer (L) in a Tyvek®¹ pouch, two separate IFUs for each respective device (PerQseal+ and PerQseal Introducer) and Quick reference guide. The PerQseal+ foil pouch contains a Tyvek® pouch which in turn contains the PerQseal+ Closure Device, housed in a thermoformed tray. Note the foil pouch is a non-sterile moisture barrier. The PerQseal Introducer pouch contains the Introducer-sheath and Dilator mounted on a support card. **The Tyvek® pouches provide the sterile barrier for their respective contents.** Contents are non-pyrogenic and supplied STERILE. All packaging not contained within the Tyvek® pouches is non-sterile. The PerQseal+ Closure Device is sterilised using an ethylene oxide (EO) process.

Please refer to the end of this booklet for graphical symbols for medical device labelling.

Store in a cool, dry, dark place.

Vivasure Medical Limited
 Parkmore Business Park West
 Galway
 Ireland, H91 V3KP
 www.vivasuremedical.com

Intended Purpose

The PerQseal+ device is intended for percutaneous vascular closure following interventional therapeutic endovascular procedures.

Indications for use

The PerQseal+ Closure Device shall be indicated for the percutaneous sealing of a common femoral arteriotomy in patients following interventional therapeutic endovascular procedures created with sheath sizes 14 - 22 F (arteriotomy up to 26 F).

The PerQseal+ Closure Device is indicated for use with 0.014" (0.356 mm) guidewires of minimum length 150 cm only.

The PerQseal+ Closure Device is indicated for use with the 'L' PerQseal Introducer.

¹ Tyvek® is a trademark of DuPont.

Clinical Benefits

- Less invasive percutaneous sealing of arteriotomy compared to surgical cut-down
- Implant is fully bioabsorbed, leaving nothing permanent behind
- Minimal pain and discomfort associated with use of VCD
- Minimisation of secondary interventions to control haemostasis
- Minimal procedural steps required to achieve haemostasis
- Percutaneous closure leads to shorter overall procedure time
- VCD has lower major complications rates than alternative therapy of cut-down and sutured closure of femoral arterial puncture sites in the indicated range
- Safe and Effective sealing of the puncture site for subjects treated with anticoagulation therapy, antiplatelet agents, intravenous glycoprotein IIb /IIIa inhibitors, or thrombolytic agents
- Delivered and deployed at the conclusion of the primary procedure (minimisation of steps for access in emergency procedures)
- Arterial wire access maintained throughout the device delivery
- Minimisation of the temporary disruption of arterial flow, which occurs with arterial clamping during surgical closure
- Reduction in scarring compared to surgical cut-down (as VCD use is percutaneous)

The Frontier V clinical study was performed to confirm safety and performance of the PerQseal+ Closure Device and PerQseal L Introducers to percutaneously close femoral artery punctures and to induce arterial haemostasis in patients undergoing endovascular procedures requiring an arteriotomy up to 26 F. The Frontier V study was a prospective, multicentre, non-randomised study to investigate the safety and performance of the PerQseal+. The study was not blinded prior to, during or post the procedure. Patients were not excluded on the basis of age, race, concomitant therapy or co-existing diseases. Primary Endpoint (Safety): Incidence of major vascular access site complications related to the PerQseal+ up to 1-month from implantation, is no worse than the performance goal for safety estimated from literature (in an equivalent patient population) for percutaneous vascular closure. Performance: Assessed by technical success rate for the PerQseal+ at discharge is no worse than success rate estimated from literature for alternative large hole closure devices from an equivalent patient population.

Study Population (4 EU sites)	Subjects N (Closures)	Male, N (%)	Age, year
	53 (61)	34 (64.2)	78.7 (±7.8)
Primary Procedure	TAVR	EVAR	
	38	15	
Follow Up	1 month	3 months	
	46	32	
Results	PerQseal+ Related Major Complication Rate	Technical Success Rate	
	0%	97%	

Conclusion: Results from the Frontier V clinical study have demonstrated the safety and efficacy of the PerQseal+ closure device for the percutaneous closure of large-bore femoral arteriotomies.

Some of the limitations of the Frontier V clinical study include:

- The patient population was selected based on the inclusion/exclusion criteria. These exclusions included a restriction of the Sheath to Femoral Artery Ratio greater than 1.05 and anterior or circumferential calcification at or near the access site.
- There are a relatively small number of subjects with a BMI greater than 30 in the study population
- There was a relatively small number of subjects enrolled in the study with a primary sheath size greater than 19 F

Post Market Clinical Follow-up studies will continue in order to confirm the safety and efficacy of the PerQseal+ in a larger patient population

The Summary of Safety and Clinical Performance of PerQseal+ closure device is available in the European database on medical devices (Eudamed), refer to <https://ec.europa.eu/tools/eudamed> and search by Basic UDI-DI (provided on page 2).

Contraindications

There are no contraindications to the use of this device. Attention is drawn to the warnings, precautions, and special patient populations. Furthermore, when used in accordance with these instructions there are no known interactions between the implant material and/or its products of metabolism and other medical devices, medicinal products and/or other substances expected to be used in conjunction with the PerQseal+ closure device.

Warnings

- Contents supplied STERILE. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilise. Reuse, reprocessing or resterilisation may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
- Do not use the Product if the packaging has been previously opened or damaged, or if any of the components appear to be damaged or defective. This may result in patient injury, illness or death.
- Discard product if 'Use By date' has passed or 'Use By date' is illegible or labelling is incomplete.
- Do not use the PerQseal+ closure device if the puncture site has been made above the most inferior border of the inferior epigastric artery (IEA) and/or above the inguinal ligament based upon bony landmarks, as this may result in retroperitoneal bleeding or haematoma. Perform a femoral angiogram or ultrasound to verify location of the puncture site.
- Do not use the PerQseal+ closure device if the puncture site has been made in the profunda femoris artery (PFA), superficial femoral artery (SFA) or located less than 15 mm proximal to the bifurcation of the profunda femoris and superficial femoral artery, since such puncture sites may result in incorrect implant position resulting in a pseudoaneurysm, intimal dissection, or vessel occlusion. Perform an oblique femoral angiogram or ultrasound to identify the bifurcation and verify location of the puncture site is within the common femoral artery.
- Do not use the PerQseal+ closure system if the percutaneous tissue tract is expected to be greater than 10 cm.

Precautions

- This product is intended for use by physicians experienced in therapeutic endovascular techniques.
- Do not advance, torque or withdraw the PerQseal+ closure device against excessive resistance. Excessive force used to advance, torque or withdraw the PerQseal+ closure device should be avoided, as this may lead to significant vessel damage and/or breakage of the device.
- Remove the PerQseal+ closure device from the patient if it is difficult to reposition to the arteriotomy.
- Remove the PerQseal+ closure device from the patient if it is difficult to deploy.
- Do not use the PerQseal+ closure device where bacterial contamination of the procedural sheath or surrounding tissues may have occurred as this may result in an infection.
- Do not use the PerQseal+ closure device if the puncture is through the posterior wall or if there are multiple punctures, since such punctures may result in a retroperitoneal bleeding.
- Do not use the PerQseal+ closure device if there is a vascular graft or stent within 20 mm of the access site.
- Do not modify the PerQseal+ closure device.
- As with all percutaneous based procedures, infection is a possibility. Observe sterile techniques at all times when using the PerQseal+ closure device. Employ appropriate groin management post procedure and post hospital discharge to prevent infection.

- If a patient has had a procedural sheath left in place for longer than 8 hours, consider the use of prophylactic antibiotics before insertion of the PerQseal+ closure device.
- Do not use the PerQseal+ closure device subsequent to or in conjunction (in the ipsilateral limb) with any other closure device or techniques, e.g. suturing.
- Following implantation of the PerQseal+ closure device, avoid vascular access at the same location within 90 days.
- Consideration should be given to the access site: tortuosity; disease state (including calcification, plaque and thrombus); and arterial diameter relative to the largest sheath size outer-diameter to ensure a successful closure and to minimise vascular complications.
- Excessive compression may damage the implant and inadvertently cause bleeding.

Special Patient Populations

The safety of the PerQseal+ closure device has not been established in the following patient populations:

- Patients with evidence of systemic bacterial or cutaneous infection, including groin infection.
- Patients suffering with known bleeding diathesis (including severe liver disease), definite or potential coagulopathy, platelet count < 100,000/µl or patients on long term anticoagulants with an INR greater than 2 at time of procedure or known type II heparin-induced thrombocytopenia.
- Patients with evidence of previous groin surgery within the region of the ipsilateral access.
- Patients with severe; claudication or peripheral vascular disease (e.g. Rutherford category 3 or greater or ABI < 0.5), documented untreated iliac artery diameter stenosis > 50% or previous bypass surgery/stent placement in the common femoral artery of ipsilateral limb.
- Patients with known allergy to any of the materials used in the PerQseal+ or PerQseal Introducer.
- Patients who have undergone a percutaneous procedure using a non-absorbable vascular closure device (excluding suture mediated) for haemostasis in the ipsilateral target leg.
- Patients with undergone a percutaneous procedure in the ipsilateral leg, within the previous 30 days.
- Patients with undergone a percutaneous procedure using an absorbable intravascular closure device for haemostasis, in the ipsilateral leg, within the previous 90 days.
- Patients with evidence of arterial diameter stenosis > 20% or anterior or circumferential calcification within 20 mm proximal or distal to target arteriotomy site.
- Females who were pregnant or lactating or not taking adequate contraceptives.
- Patients with a lower extremity amputation from the ipsilateral or contralateral limb.
- Patients with significant blood loss/transfusion (defined as requiring transfusion of 4 or more units of blood products) during index procedure or within 30 days prior to index procedure.
- Patients with unstable blood pressure or heart rate, symptoms of shock or systolic blood pressure greater than 160 mmHg at the time of closure.
- Patients with an activated clotting time (ACT) > 350 seconds immediately prior to sheath removal or patients where the ACT measurements are expected to be > 350 seconds for more than 24 hours after closure procedure.
- Patients with an acute haematoma of any size, arteriovenous fistula or pseudoaneurysm at the target access site; or any evidence of arterial laceration or dissection within the external iliac or femoral artery before the use of the PerQseal+ Closure Device.
- Patients on thrombolytic agents within 24 hours prior to or during the endovascular procedure which caused their fibrinogen level to be < 100 mg/dl.
- Patients under the age of 18 years.

Procedure

The PerQseal+ closure device is to be used only by a licensed physician possessing documented instruction and training by Vivasure Medical in the use of the PerQseal+ closure device. Exercise care during PerQseal+ closure device handling to reduce the possibility of accidental device breakage.

It is recommended to perform angiography in both the Anterior Posterior projection and Oblique projection (Right Anterior Oblique – 30 degrees from the right limb, Left Anterior Oblique – 30 degrees from the left limb) prior to ipsilateral access to identify the bifurcation and ensure that the location of the puncture will be in the common femoral artery, at the 12 o'clock position.

It is recommended that ultrasound or angiographic guided puncture / micropuncture is used for access via the ipsilateral limb, to ensure a centre stick, away from anterior calcification, at least 15 mm proximal to the SFA/PFA bifurcation and access is within the common femoral artery. Doppler ultrasound may be used at the completion of Procedural STEP 2 to ensure implant is correctly positioned juxtaposed to the arteriotomy prior to completing device deployment.

The PerQseal+ closure device should be used within four hours of opening the foil pouch.

If the PerQseal+ closure device cannot be positioned juxtaposed to the arteriotomy, due to improper orientation or patient vascular anatomy, the implant and delivery system should be withdrawn from the patient. Haemostasis can then be achieved by employing good clinical judgement.

Examination of Products

1. Carefully inspect the packaging to ensure there is no damage.
2. Check all Labelling to ensure the product is within its use by date.

Opening the packaging

1. Remove the foil pouch from the procedure pack.
2. Inspect the foil pouch to ensure that it has not been breached.
3. Retain the implant card for provision to the patient post procedure.
4. Tear open the top of the foil pouch from the side at the notched position.
5. Open the end of the foil pouch and remove the Tyvek® pouch.
6. Inspect the Tyvek® pouch to ensure the sterile barrier has not been breached.
7. While observing aseptic techniques, open the Tyvek® pouch by peeling apart the end with the chevron seal (the end with the 'v-shape' seal).
8. The device must be used within 4 hours of opening the foil pouch.
9. Present the exposed tray for aseptic transfer into the sterile field.
10. Have a person in the sterile field remove the tray containing the PerQseal+ closure device from the Tyvek® pouch.

Preparation

1. Insert and position the PerQseal L Introducer-sheath as per PerQseal Introducer IFU.
2. Using standard clinical practice insert a 0.014" (0.356 mm) guidewire with at least 50 cm of available guidewire protruding from the Introducer-sheath.
Note: If a vascular graft is present in the expected 0.014" guidewire path consider delivering the 0.014" guidewire via a pigtail catheter once the dilator has been removed from the Introducer-sheath.
3. Ensure that the Introducer-sheath is positioned at least 4 graduations into the vessel as per PerQseal Introducer IFU. The Introducer-sheath should be at a low angle, of less than 40° (ideally <30°) relative to the longitudinal axis of the femoral artery.
4. Perform an ACT reading to ensure that ACT is less than 350 seconds.

5. It is recommended to perform an ipsilateral site angiogram to confirm the arteriotomy puncture is within the common femoral artery and to assess presence of:
- Dissection
 - Stenosis (greater than 20% within 20 mm of the arteriotomy)
 - Perforation other than the access arteriotomy

Note: If any of the above is confirmed by the angiogram, do not use the PerQseal+ closure device.

Tray Loading

- A. Whilst the tray is lying flat on the sterile preparation area, fill the recess at the end of the tray (indicated on the tray with “A”) with sterile non-heparinised saline solution (approximately 20 cc) to submerge the implant, Figure 2.


Figure 2. Fill tray recess with non-heparinised saline.

- B. Pull off the small Removable tray lid, indicated with a “B” on the tray, by lifting the tab as shown in Figure 3 until completely removed.

Note: do not remove any other lids.


Figure 3. Remove Removable Tray Lid.

- C. Holding the tray steady, grip the rear of the handle (release cam) and pull until a firm stop is reached (indicated on the tray with “C”), Figure 4. The implant is now folded in the loading cannula tube.


Figure 4. Pull device handle back to firm stop.

- D. Flip the funnel lever to a vertical position to release the delivery device from the loading funnel (indicated on the tray with “D”), Figure 5.


Figure 5. Flip the funnel lever to release device from tray.

- Continue to pull back and withdraw the PerQseal+ closure device from the tray.
- Visually check the integrity of the implant (within the loading cannula); there should be no creases, folds or tears, see Figure 6.


Figure 6. Inspection of Implant (within loading cannula).

- The PerQseal+ closure device is now ready for use.

Note: Use the PerQseal+ closure device within 30 minutes of completing this step.

Procedural Steps:

Insert the Device

Load Guidewire:

- Use non-heparinised saline to lubricate the guidewire during insertion. Insert the guidewire into the implant within the loading cannula, Figure 7 continue to insert the guidewire until it emerges from the rear of the PerQseal+ closure device handle. Hold the PerQseal Introducer-sheath at a low angle relative to the longitudinal axis of the artery (e.g. 20°) and keep the Introducer-sheath as straight as practical, to avoid kinking.


Figure 7. Loading guidewire into device cannula.

First Lock:

- Advance the Loading cannula over the guidewire towards the Introducer-sheath hub, Figure 8. If difficulty is experienced tracking over the guidewire, ensure the guidewire is not kinked and that it is free of coagulated blood and thrombus.

Note 1: Hold the guidewire at the rear of the handle during insertion to ensure the PerQseal+ closure device travels over the wire and do not inadvertently insert additional wire into the patient during PerQseal+ closure device insertion/advancement.

Note 2: Only grip the back portion of the Loading cannula during insertion as indicated in Figure 7.


Figure 8. Advance the PerQseal+ closure device to the Introducer-sheath.

Maintaining the position of the Introducer-sheath insert the loading cannula fully into the Introducer-sheath hub. The loading cannula cap should be fully seated and locked into the Introducer-sheath hub, (with a click). The green marking being aligned confirms this step is successfully completed, Figure 9.


Figure 9. Seating Loading cannula into Introducer-sheath

Second Lock:

- Holding the Introducer-sheath hub in one hand and the device handle in the other hand, carefully advance the device handle. Ensure the indicators on the Introducer-sheath and handle are in alignment and the loading cannula/sheath hub assembly locks (with a click) into the device handle. Note the green indicator in the window. This confirms that the PerQseal+ closure device handle has been successfully locked. Figure 10.

Note: During insertion, maintain the guidewire stationary at the rear of the handle. Ensure the position of the sheath in the artery remains unchanged. At this stage the handle, loading cannula and sheath are all locked together as one.


Figure 10. Insertion of PerQseal+ closure device into Introducer-sheath.

Check the orientation/alignment of the PerQseal+ closure device by ensuring the handle and Introducer-sheath markings are facing vertically upwards away from the patient (see Figure 11). If required, rotate the device (using the handle) to achieve the correct alignment.


Figure 11. Correct orientation of PerQseal+ closure device prior to deployment.

Grip the distal end of the PerQseal+ closure device handle (closest to patient). Stabilise the device in position by supporting this hand on the patient, Figure 12.

STEP 1 Unsheathing the implant

- Place hand on the sheath cam, marked '1' and rotate this grip section in the indicated direction 180° to align the dots. Ensure the hand gripping the distal end of the handle remains fixed and the handle remains stationary during this rotation, Figure 12. The implant is now exposed within the artery.

Note 1: the Introducer-sheath will withdraw from the patient into the handle of the PerQseal+ closure device during this step to expose the implant in the artery. It is important that the handle does not advance toward the patient during rotation of the cam.

Note 2: If high resistance to rotation is felt prior to reaching the 180° rotation relieve the rotation pressure and try again. Do not rotate the cam backward.

Note 3: Do not rotate beyond 180°, this is indicated by the alignment of green dots.


Figure 12. Step 1 of PerQseal+ closure device Actuation (Unsheathing the implant).

STEP 2 Tamponade

- While holding the device handle at a low angle (≤ 30 degrees), apply light extracorporeal pressure over the skin insertion site with the other hand. Then slowly withdraw the device from the artery, whilst watching the graduations on the Introducer-sheath, when approximately two graduations have been withdrawn, the Introducer-sheath tip will have exited the arteriotomy. Expect a reduction in tension needed to withdraw the PerQseal+ closure device, after the Introducer-sheath tip exits the arteriotomy.

Note: The digital pressure should ONLY be applied by the device operator to ensure tactile feedback occurs.

- Continue to withdraw the device a further 1.5 graduations (= 3.5 graduations total). At this point, gentle tactile resistance should be felt as a result of the implant anchoring itself against the inside lumen of the artery wall. This is the tamponade position, and the implant is now controlling the bleeding from the arteriotomy. Figure 13.

Note 1: It is important to confirm the location by gentle resistance of the implant against the arterial wall to avoid the implant being deployed within the tissue tract.

Note 2: Avoid excessive tension on the PerQseal+ closure device. Excessive tension on the PerQseal+ closure device will reduce the implant's effectiveness to achieve haemostasis or could result in the implant being pulled out of the arteriotomy.

Note 3: The PerQseal+ closure device delivery system location at the arteriotomy may be confirmed by fluoroscopy/angiography or Doppler Ultrasonography.


Figure 13. Tamponade position, implant positioned juxtaposition to arteriotomy.

- c. Slight adjustment of PerQseal+ closure device position may be necessary to achieve tamponade. Ensure that an acceptable tamponade is achieved before proceeding with deployment.

NOTE: Once an acceptable tamponade is achieved, avoid implant/sheath movement until the device is fully deployed.

STEP 3 Remove guidewire

- a. Locate indicator marked '3' and remove the guidewire completely from the patient and PerQseal+ closure device. Guidewire removal may be done by an assistant. Figure 14.


Figure 14. Removal of guidewire.

Note: The position of the PerQseal+ closure device handle should remain fixed at all times. (Such that the implant is positioned at the arteriotomy with minimal tension.)

STEP 4. Release cam lock

- a. Locate the cam lock (marked '4') and fully depress this green button. Figure 15. Ensure it is fully depressed before proceeding to step 5.


Figure 15. Release cam lock.

STEP 5 Implant Deployment and Release

- a. Place the free hand on the Release cam, marked '5' and partially rotate this grip section 90° in the direction indicated as far as the first alignment marks. Figure 16. The implant is now in the fully deployed configuration, but the implant is still attached to the device handle.

Note: Ensure that tamponade is still acceptable before proceeding with full deployment. Slight adjustment of PerQseal+ closure device position may be necessary to achieve tamponade.

- b. Fully rotate this grip section (marked '5') a further 90° in the direction indicated as far as the second alignment marks. This releases the implant from the delivery system.

Note: Ensure the hand gripping the distal end of the handle remains fixed and the handle remains stationary during both stages of Step 5.


Figure 16. Release implant.

- c. Gently withdraw the delivery system from the tissue tract. Figure 17.


Figure 17. Removal of deployed device.

- d. If suitable access is available an angiogram or duplex ultrasound may be performed immediately after PerQseal+ closure device deployment to confirm arteriotomy sealing.
Note: If haemostasis has not occurred, gentle manual compression may be applied proximal to the arteriotomy sufficient to control bleeding. Consult troubleshooting section.
- e. Use standard clinical procedures to treat any skin incisions.

Post procedural Management

Provide patient with Patient Implant Card containing all of the required information.

Before considering patient discharge, assess for the following:

- Bleeding from the access site
- Anticoagulation thrombolytic or antiplatelet therapy
- Haematoma at the access site
- Hypotension
- Any condition requiring clinical observation
- Pain, discomfort or transitory local irritation and inflammation at the puncture site.
- Local and distal pulse deficits, ischaemia.

TROUBLE SHOOTING GUIDE

- 1. Cannot insert the loading cannula into the Introducer-sheath hub:**
Ensure that the supplied PerQseal Introducer-sheath is being used and the loading cannula is not damaged.
- 2. Difficult to advance the PerQseal+ closure device into the Introducer-sheath:**
Ensure that the loading cannula is fully engaged with the Introducer-sheath hub, if not discard the PerQseal+ closure device and replace with new PerQseal+ closure device.
- 3. PerQseal+ closure device implant withdrawn from artery during positioning:**
Control bleeding by application of manual compression or cross-over balloon if *in situ*. Withdraw the PerQseal+ closure device and Introducer-sheath from the patient, as per Section 1 below (Recommended Recovery Options), ensuring to leave the guidewire *in situ*. Discard the PerQseal+ closure device and Introducer-sheath and replace with new PerQseal+ closure device.
- 4. Cannot achieve acceptable tamponade during implant positioning prior to deployment:**
 - Avoid excessive tension on the PerQseal+ closure device as excessive tension will reduce the implant's effectiveness to achieve haemostasis and control bleeding.
 - Change the angle and/or orientation of the PerQseal+ closure device (relative to the longitudinal axis of the artery) by using very small adjustments (approximately 5°).
 - If acceptable tamponade cannot be achieved do not deploy the PerQseal+ closure device. Follow the recover suggestion outlined in Section 1 below (Recommended Recovery Options).
- 5. Difficult to withdraw the guidewire from the PerQseal+ closure device:**
 - Do not force.
 - Ensure the guidewire is not snagged endovascularly.
 - Lower the angle of the PerQseal+ closure device (relative to the longitudinal axis of the artery).
 - If it is not possible to withdraw the guidewire do not deploy the PerQseal+ closure device and follow the steps outlined in Section 1 Recommended Recovery Options below.
- 6. Haemostasis not achieved post deployment of the PerQseal+ closure device:**
 - Light manual compression may be applied proximal to the implant to control bleeding until haemostasis is achieved.
 - Note: It is not recommended to apply manual pressure directly onto the location where implant is situated. Excessive compression may damage the implant and adversely cause bleeding.
 - Ensure ACT level is below 350 seconds, use clinical judgement and reduce the effects of heparin with agents such as protamine sulphate.
 - Coagulating extraarterial agents may be used.
 - Temporary use of a cross-over balloon inflated directly beneath the implant may be used.
 - Keep patient under observation until stable.

7. Difficult to withdraw the device post deployment:

The device tip may be caught on subcuticular soft tissue. Rotate the device in either direction and try to withdraw.

8. Guidewire in situ post deployment:

If the guidewire is still *in situ* post deployment, follow the recovery steps outlined in Section 4 Recommended Recovery Options.

Recommended Recovery Options

1. If it becomes necessary to remove the implant from the artery whilst still attached to the delivery system then follow these recommended recovery steps:
 - i) Control bleeding by manual compression or via cross-over balloon inflation.
 - ii) Ensure Introducer-sheath is withdrawn from the artery (bring device to **Step 2 Achieve Tamponade** stage).
 - iii) Advance the PerQseal+ closure device into the artery approximately 1 cm.
 - iv) Rotate the handle 180° so that the handle is now facing towards the patient in a downward direction.
 - v) Withdraw the handle to remove the PerQseal+ closure device implant from the artery whilst maintaining the PerQseal+ closure device at an approximate angle of 20° relative to the patient. Ensure the guidewire remains *in situ* during withdrawal.
 - vi) Re-cannulate the artery over the 0.014" guidewire with 4-9 F dilator.
 - vii) Exchange the 0.014" guidewire for a 0.035" guidewire.
 - viii) Open a new PerQseal+ closure device procedure pack.
 - ix) Remove the 4-9 F dilator and reinsert a new Introducer-sheath as per PerQseal Introducer IFU.
 - x) Exchange the 0.035" guidewire for a suitable 0.014" guidewire, ensuring approximately 50 cm length of guidewire is available from the rear of the Introducer-sheath.
 - xi) Follow the PerQseal+ closure device deployment steps outlined above.
2. If the implant becomes detached from the delivery system within the artery whilst the Introducer-sheath and guidewire are *in situ* use the following recovery steps to retrieve the implant:
 - i) Withdraw the delivery system and the Introducer-sheath and ensure the guidewire remains *in situ*.
 - ii) Re-cannulate the artery over the 0.014" guidewire with a 4-9 F dilator.
 - iii) Insert a 0.035" guidewire past the implant.
 - iv) Open a new PerQseal+ device Procedure Pack and insert Introducer as per PerQseal Introducer IFU.
 - v) Remove the dilator.
 - vi) Insert a snare beyond the implant to snare the original 0.014" guidewire.
 - vii) Withdraw the snare, original 0.014" guidewire and trapped implant into the Introducer-sheath.
 - viii) Ensure the implant is within the Introducer-sheath and remove the snare, original guidewire, trapped implant and Introducer-sheath from the patient, whilst maintaining the position of the second 0.035" guidewire *in situ*.
 - ix) Insert PerQseal Introducer and deploy the new PerQseal+ closure device as per instructions above.
3. For recovery of the bioabsorbable implant from within the artery post Introducer-sheath and guidewire removal it is recommended to locate the implant via angiography. Based on the implant's location and the level of perfusion downstream of the implant, use clinical judgement to decide an appropriate strategy for recovery.
4. If the device has been fully deployed whilst the 0.014" guidewire is still *in situ* use the following recovery steps to retrieve the guidewire.
 - i) Remove the sheath and delivery system from the guidewire.
 - ii) Very slowly increase tension on the guidewire, at the same angle as deployment or approx. 30° to the artery, until movement of the guidewire occurs. Simultaneously apply light extracorporeal pressure slightly distal to the position of the implant.
 - iii) Continue to withdraw the guidewire whilst maintaining light extracorporeal pressure until fully removed.
 - iv) Inspect the guidewire to ensure the entire guidewire has been removed from the patient.
 - v) Expect some bleeding from the guidewire port in the implant. This can be managed by applying sufficient manual compression to control bleeding. Note: excessive compressing will damage the implant and adversely cause bleeding.
 - vi) Keep the patient under observation until stable.

Adverse Events & Risks

The following are anticipated complications associated with use of the PerQseal+ closure device:

- Bruising, oozing at the puncture site
- Transient extravasation
- Late bleed
- Retroperitoneal bleeding
- Haematoma
- Pain, discomfort or transitory local irritation and inflammation at the puncture site
- A mild local inflammatory reaction to the implant material (related to the ongoing biodegradation process) within the first 90 days. The breakdown products of the implant material are glyoxylate or glycine which is extracted via urine or as carbon dioxide respectively
- Access-site related nerve injury, vascular spasm
- Local and/or distal pulse deficits, ischemia
- Stenosis at the closure site
- Infection
- Arterial or deep vein thrombosis
- Pseudoaneurysm
- Arteriovenous fistula
- Vasovagal response

Other events that could possibly occur include:

- Occlusive intraluminal thrombus and/or emboli formation at the implant site
- Sinus formation
- Compartment syndrome
- Iliac or femoral artery dissection
- Partial or complete occlusion of the artery or peripheral arteries
- Embolisation of the intraarterial implant components
- Damage of the arterial wall (perforation or rupture)
- Allergic or adverse foreign body reaction
- Corrective intervention due to any of the above complications such as surgery and/or transfusion

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the procedure was performed.

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