PerQseal[®] Elite

CLOSURE DEVICE

Notice: Use of the PerQseal[®] Elite Closure Device should be performed only by physicians who have successfully completed Vivasure approved training

En: INSTRUCTIONS FOR USE	2
De: GEBRAUCHSANWEISUNG	14
Fr : MODE D'EMPLOI	28
NL: GEBRUIKSAANWIJZING	42
Da: BRUGSANVISNING	56
Es: INSTRUCCIONES DE USO	70
El: ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	84

PerQseal[®] Elite

INSTRUCTIONS FOR USE

To ensure proper deployment and use of this device and to prevent injury to patients, read all information contained in these instructions

For assistance, contact:
Vivasure Medical customer Service:

Telephone Number (9:00 – 17:00 GMT Mon to Fri)	+353 91 395 440
Basic UDI-DI	539153318DP4FA1VCDV7

Warranty

VIVASURE warrants that reasonable care has been used in the design and manufacture of this device. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilisation of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond VIVASURE's control directly affect the device and the results obtained from its use. VIVASURE shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. VIVASURE neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. **VIVASURE assumes no liability with respect to devices reused, reprocessed or resterilised and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such devices.**

TABLE OF CONTENTS

ENGLISH/ENGLISH	
WARRANTY	2
DEVICE DESCRIPTION	3
HOW SUPPLIED	4
INTENDED PURPOSE	4
INDICATIONS FOR USE	4
CLINICAL BENEFITS	4
CONTRAINDICATIONS	5
WARNINGS	5
PRECAUTIONS	5
SPECIAL PATIENT POPULATIONS	5
PROCEDURE	6
POST PROCEDURAL MANAGEMENT	11
TROUBLE SHOOTING GUIDE	12
RECOMMENDED RECOVERY OPTIONS	12
ADVERSE EVENTS & RISKS	13
GRAPHICAL SYMBOLS FOR MEDICAL DEVICE LABELLING	111

Device Description

The PerQseal Elite closure device is a Vascular Closure Device (VCD), designed to percutaneously deliver a bioabsorbable vascular implant (with the presence of an abluminal surface matrix) to seal femoral artery puncture sites following endovascular procedures.

The implant has both intraarterial and extraarterial components. The implant is inserted into the artery through an Introducer-sheath and over a guidewire. The implant is designed to be positioned relative to the arteriotomy and effect a seal. The implant has passed design validation testing in pre-clinical and clinical studies in the presence of procedural fluids used in the handling and preparation of the implant without evidence of interaction with the absorbable material.

The PerQseal Elite closure device has a delivery system for the percutaneous delivery of the implant to the arteriotomy. The delivery system has a handle at its proximal end and a flexible shaft, which attaches to the implant at the distal end.

Refer to Figure 1a showing the delivery system and implant, indicating the key features of the PerQseal Elite closure device. The delivery system is designed to ensure single use of the device.

Refer to Figure 1b showing enlarged image of the implant in deployed configuration and identifying each component.


Figure 1a. PerQseal[®] Elite Closure Device

1. Handle	6. Release cam (Back cam)
2. Shaft	7. Loading cannula
3. Implant	8. Fixation cam
4. Sheath cam	9. Funnel lever
5. Cam lock	10. Removable tray lid


Figure 1b. Implant (Deployed configuration)

1. Patch	3. Scaffold
2. External Fixation	4. Guidewire Pin

The PerQseal Elite Closure Device will be used in conjunction with the following:

- PerQseal[®] Elite Introducer (Introducer-sheath and Dilator)
- A 0.035" (0.89mm) guidewire of minimum length 150 cm – not provided
- 20 ml syringe – not provided

Note: The term 'Device' is defined as the entire VCD system including delivery system, and implant. The term 'Product' is defined as the entire VCD system including delivery system, implant, Instructions for Use (IFU), labelling and packaging.

Implant Composition: < 223 mg Polydioxanone. Polydioxanone is a synthetic, absorbable, low reactivity polymer, degradation of which occurs along the polymer back bone by hydrolysis of the ester bonds. This occurs randomly along the polymer chain cutting the bonds and converting the long polymer chains into (shorter fragments) of natural monomeric acids found in the body. The polymer is thus degraded to non-toxic, low molecular weight residues capable of being eliminated from the body by normal metabolic pathways. The expected lifetime of the device is ~30 days at which point the implant begins to lose mechanical functionality. Total absorption of the implant is expected at ~ 180-days.

Maximum of one implant per patient.

This implant has been evaluated and determined to be MR Safe. Patients implanted with this device may safely undergo Magnetic Resonance Imaging (MRI) procedures.

Delivery System Composition: The delivery system materials the patient could be exposed to include stainless steel (contains Nickel), nylon, acrylonitrile butadiene styrene (ABS), polyether ether ketone (PEEK).

How Supplied

The packaged device is supplied in the 'PerQseal Elite' procedure pack. Each procedure pack contains the 'PerQseal Elite Closure Device' in a foil pouch, the 'PerQseal Elite Introducer' in a Tyvek^{®1} pouch, two separate IFUs for each respective device (PerQseal Elite Closure Device and PerQseal Elite Introducer) a Patient's Implant Card and Quick reference guide. The foil pouch contains a Tyvek pouch which in turn contains the PerQseal Elite Closure Device, housed in a thermoformed tray. Note the foil pouch is a non-sterile moisture barrier. The PerQseal Elite Introducer pouch contains the Introducer-sheath and Dilator mounted on a support card. **The Tyvek pouches provide the sterile barrier system for their respective contents.** Contents are non-pyrogenic and supplied STERILE. All packaging not contained within the Tyvek pouches is non-sterile. The PerQseal Elite Closure Device is sterilised using an ethylene oxide (EO) process.

Please refer to the end of this booklet for graphical symbols for medical device labelling.

Store in a cool, dry, dark place.

Vivasure Medical Limited
Parkmore Business Park West
Galway
Ireland
H91 V3KP.
www.vivasuremedical.com

Intended Purpose

The PerQseal Elite Vascular Closure Device is intended for percutaneous common femoral arteriotomy closure following interventional therapeutic endovascular procedures.

Indications for use

The PerQseal Elite Closure Device shall be indicated for the percutaneous sealing of a common femoral arteriotomy in patients following interventional therapeutic endovascular procedures created with sheath sizes 14 - 22 F (arteriotomy up to 26 F).

The PerQseal Elite Closure Device is indicated for use with 0.035" (0.89 mm) guidewires of minimum length 150 cm only.

The PerQseal Elite Closure Device shall be indicated for use with the PerQseal Elite Introducer.

Intended Patient Population

- Individuals over 18 years of age.
- Females who are not pregnant or lactating.
- Patients that are not vulnerable.

Intended User

The intended users of the PerQseal Elite include, but are not limited to, vascular surgeons, interventional radiologists, interventional cardiologists, cardiothoracic surgeons and angiologists. All users are expected to be medically trained and qualified practitioners.

Clinical Benefits

- Less invasive percutaneous sealing of arteriotomy compared to surgical cut-down.
- Implant is fully bioabsorbed within 180 days, leaving nothing permanent behind
- Minimal pain and discomfort associated with use of VCD
- Minimal procedural steps required to achieve haemostasis
- Percutaneous closure leads to shorter overall procedure time
- VCD has lower major complications rates than alternative therapy of cut-down and sutured closure of femoral arterial puncture sites in the indicated range
- Safe and Effective sealing of the puncture site for subjects treated with anticoagulation therapy, antiplatelet agents, intravenous glycoprotein IIb /IIIa inhibitors, or thrombolytic agents
- Reduction in scarring compared to surgical cut-down (as VCD use is percutaneous)

The Frontier V clinical study was performed to confirm safety and performance of an equivalent device, the PerQseal+ Closure Device and PerQseal L Introducers, to percutaneously close femoral artery punctures and to induce arterial haemostasis in patients undergoing endovascular procedures requiring an arteriotomy up to 26 F. The Frontier V study was a prospective, multicentre, non-randomised study to investigate the safety and performance of the PerQseal+. The study was not blinded prior to, during or post the procedure. Patients were not excluded on the basis of age, race, concomitant therapy or co-existing diseases.

Primary Endpoint (Safety): Incidence of major vascular access site complications related to the PerQseal+ up to 1-month from implantation, is no worse than the performance goal for safety estimated from literature (in an equivalent patient population) for percutaneous vascular closure.

Performance: Assessed by technical success rate for the PerQseal+ at discharge is no worse than success rate estimated from literature for alternative large hole closure devices from an equivalent patient population.

Study Population (7 EU sites)	Subjects N (Closures)	Male, N (%)	Age, years
	91 (103)	63 (69%)	78.2 (±8.1)
Primary Procedure	TAVR	EVAR	
	56	35	
Follow Up	1 month	3 months	
	87	84	
Results	PerQseal+ Related Major Complication Rate, N (%)	Technical Success Rate	Time to haemostasis (minutes) [Mean ± SD, Median]
	1 (1%)	95%	4.0 ± 5.0, 1.0

Conclusion: Results from the Frontier V clinical study have demonstrated the safety and efficacy of the PerQseal+ closure device for the percutaneous closure of large-bore femoral arteriotomies.

¹ Tyvek[®] is a trademark of DuPont.

Some of the limitations of the Frontier V clinical study include:

- The patient population was selected based on the inclusion/exclusion criteria. These exclusions included a restriction of the Sheath to Femoral Artery Ratio greater than 1.05 and anterior or circumferential calcification at or near the access site.
- There are a relatively small number of subjects with a BMI greater than 30 in the study population
- There was a relatively small number of subjects enrolled in the study with a primary sheath size greater than 19 F

Post Market Clinical Follow-up studies will continue in order to confirm the safety and efficacy of the PerQseal+ in a larger patient population, and to confirm the safety and efficacy of the PerQseal Elite Closure Device.

The Summary of Safety and Clinical Performance of PerQseal Elite closure device is available in the European database on medical devices (Eudamed), refer to <https://ec.europa.eu/tools/eudamed> and search by Basic UDI-DI (provided on page 2) or by request from Vivasure Medical Ltd. (contact details on page 2).

Contraindications

There are no contraindications to the use of this device. Attention is drawn to the warnings, precautions, and special patient populations. Furthermore, when used in accordance with these instructions there are no known interactions between the implant material and/or its products of metabolism and other medical devices, medicinal products and/or other substances expected to be used in conjunction with the PerQseal Elite closure device.

Warnings

- Contents supplied STERILE. Do not use if sterile barrier system is damaged.
- For single use only. Do not reuse, reprocess or resterilise. Reuse, reprocessing or resterilisation may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
- Do not use the Product if the packaging has been previously opened or damaged, or if any of the components appear to be damaged or defective. This may result in patient injury, illness or death.
- Discard product if 'Use By date' has passed or 'Use By date' is illegible or labelling is incomplete.
- Do not use the PerQseal Elite closure device if the puncture site has been made above the most inferior border of the inferior epigastric artery (IEA) and/or above the inguinal ligament based upon bony landmarks, as this may result in retroperitoneal bleeding or haematoma. Perform a femoral angiogram or ultrasound to verify location of the puncture site.
- Do not use the PerQseal Elite closure device if the puncture site has been made in the profunda femoris artery (PFA), superficial femoral artery (SFA) or located less than 15 mm proximal to the bifurcation of the profunda femoris and superficial femoral artery, since such puncture sites may result in incorrect implant position resulting in a pseudoaneurysm, intimal dissection, or vessel occlusion. Perform an oblique femoral angiogram or ultrasound to identify the bifurcation and verify location of the puncture site is within the common femoral artery.
- Do not use the PerQseal Elite closure system if the percutaneous tissue tract is expected to be greater than 10 cm.

Precautions

- This product is intended for use by physicians experienced in therapeutic endovascular techniques.
- Do not advance, torque or withdraw the PerQseal Elite closure device against excessive resistance. Excessive force used to advance, torque or withdraw the PerQseal Elite closure device should be avoided, as this may lead to significant vessel damage and/or breakage of the device.
- Remove the PerQseal Elite closure device from the patient if it is difficult to reposition to the arteriotomy.
- Remove the PerQseal Elite closure device from the patient if it is difficult to deploy.
- Do not use the PerQseal Elite closure device where bacterial contamination of the procedural sheath or surrounding tissues may have occurred as this may result in an infection.
- Do not use the PerQseal Elite closure device if the puncture is through the posterior wall or if there are multiple punctures, since such punctures may result in a retroperitoneal bleeding.
- Do not use the PerQseal Elite closure device if there is a vascular graft or stent within 20 mm of the access site.
- Do not modify the PerQseal Elite closure device.
- As with all percutaneous based procedures, infection is a possibility. Observe sterile techniques at all times when using the PerQseal Elite closure device. Employ appropriate groin management post procedure and post hospital discharge to prevent infection.
- If a patient has had a procedural sheath left in place for longer than 8 hours, consider the use of prophylactic antibiotics before insertion of the PerQseal Elite closure device.
- Do not use the PerQseal Elite closure device subsequent to or in conjunction (in the ipsilateral limb) with any other closure device or techniques, e.g. suturing.
- Following implantation of the PerQseal Elite closure device, avoid vascular access at the same location within 90 days.
- Consideration should be given to the access site: tortuosity; disease state (including calcification, plaque and thrombus); and arterial diameter relative to the largest sheath size outer-diameter to ensure a successful closure and to minimise vascular complications.
- Excessive compression may damage the implant and inadvertently cause bleeding.

Special Patient Populations

The safety of the PerQseal Elite closure device has not been established in the following patient populations:

- Patients with evidence of systemic bacterial or cutaneous infection, including groin infection.
- Patients suffering with known bleeding diathesis (including severe liver disease), definite or potential coagulopathy, platelet count < 100,000/ μ l or patients on long term anticoagulants with an INR greater than 2 at time of procedure or known type II heparin-induced thrombocytopenia.
- Patients with evidence of previous groin surgery within the region of the ipsilateral access.
- Patients with severe; claudication or peripheral vascular disease (e.g. Rutherford category 3 or greater or ABI < 0.5), documented untreated iliac artery diameter stenosis > 50% or previous bypass surgery/stent placement in the common femoral artery of ipsilateral limb.
- Patients with known allergy to any of the materials used in the PerQseal Elite or PerQseal Introducer.
- Patients who have undergone a percutaneous procedure using a non-absorbable vascular closure device (excluding suture mediated) for haemostasis in the ipsilateral target leg.
- Patients who have undergone a percutaneous procedure in the ipsilateral leg, within the previous 30 days.
- Patients who have undergone a percutaneous procedure using an absorbable intravascular closure device for haemostasis, in the ipsilateral leg, within the previous 90 days.
- Patients with evidence of arterial diameter stenosis > 20% or anterior or circumferential calcification within 20 mm proximal or distal to target arteriotomy site.
- Females who were pregnant or lactating or not taking adequate contraceptives.
- Patients with a lower extremity amputation from the ipsilateral or contralateral limb.
- Patients with significant blood loss/transfusion (defined as requiring transfusion of 4 or more units of blood products) during index procedure or within 30 days prior to index procedure.
- Patients with unstable blood pressure or heart rate, symptoms of shock or systolic blood pressure greater than 160 mmHg at the time of closure.
- Patients with an activated clotting time (ACT) > 350 seconds immediately prior to sheath removal or patients where the ACT measurements are expected to be > 350 seconds for more than 24 hours after closure procedure.
- Patients with an acute haematoma of any size, arteriovenous fistula or pseudoaneurysm at the target access site; or any evidence of arterial laceration or dissection within the external iliac or femoral artery before the use of the PerQseal Elite Closure Device.
- Patients on thrombolytic agents within 24 hours prior to or during the endovascular procedure which caused their fibrinogen level to be < 100 mg/dl.
- Patients under the age of 18 years.

PROCEDURE

The PerQseal Elite closure device is to be used only by a licensed physician possessing documented instruction and training by Vivasure Medical in the use of the PerQseal Elite closure device. Exercise care during PerQseal Elite closure device handling to reduce the possibility of accidental device breakage.

It is recommended to perform angiography in both the Anterior Posterior projection and Oblique projection (Right Anterior Oblique – 30 degrees from the right limb, Left Anterior Oblique – 30 degrees from the left limb) prior to ipsilateral access to identify the bifurcation and ensure that the location of the puncture will be in the common femoral artery, at the 12 o'clock position.

It is recommended that ultrasound or angiographic guided puncture / micropuncture is used for access via the ipsilateral limb, to ensure a centre stick, away from anterior calcification, at least 15 mm proximal to the SFA/PFA bifurcation and access is within the common femoral artery. Doppler ultrasound may be used at the completion of Procedural STEP 2 to ensure implant is correctly positioned juxtaposed to the arteriotomy prior to completing device deployment.

The PerQseal Elite closure device should be used within four hours of opening the foil pouch.

If the PerQseal Elite closure device cannot be positioned juxtaposed to the arteriotomy, due to improper orientation or patient vascular anatomy, the implant and delivery system should be withdrawn from the patient. Haemostasis can then be achieved by employing good clinical judgement.

Ensure there is no guidewire or sheath/catheter crossing the arteriotomy or within 20 mm distal or proximal to the access site prior to PerQseal Elite use.

Examination of Packaging & Labelling

1. Inspect visually the sterile barrier system (Tyvek Pouch) to ensure there is no damage or breaches of packaging integrity.
2. Check all Labelling to ensure the product is within its use by date.

Opening the packaging

1. Remove the foil pouch from the procedure pack.
2. Inspect the foil pouch to ensure that it has not been breached.
3. Retain the implant card for provision to the patient post procedure.
4. Tear open the top of the foil pouch from the side at the notched position.
5. Open the end of the foil pouch and remove the Tyvek pouch.
6. Inspect visually the sterile barrier system (Tyvek Pouch) to ensure no breaches of packaging integrity.
7. While observing aseptic techniques, open the Tyvek pouch by peeling apart the end with the chevron seal (the end with the 'v-shape' seal).
8. The device must be used within 4 hours of opening the foil pouch.
9. Present the exposed tray for aseptic transfer into the sterile field.
10. Have a person in the sterile field remove the tray containing the PerQseal Elite closure device from the Tyvek pouch.

Preparation

1. Insert and position the PerQseal Elite Introducer-sheath over the *in-situ* 0.035" guidewire as per PerQseal Elite Introducer IFU. Ensure at least 50 cm of available guidewire protruding from the Introducer-sheath.
2. Ensure that the Introducer-sheath is positioned at least 4 graduations into the vessel as per PerQseal Elite Introducer IFU. The Introducer-sheath should be at a low angle, of less than 40° relative to the longitudinal axis of the femoral artery (ideally < 30°).
3. Perform an ACT reading to ensure that ACT is less than 350 seconds.
4. It is recommended to perform an ipsilateral site angiogram to confirm the arteriotomy puncture is within the common femoral artery and to assess presence of:
 - Dissection
 - Stenosis (greater than 20% within 20 mm of the arteriotomy)
 - Perforation other than the access arteriotomy

Note: If any of the above is confirmed by the angiogram, do not use the PerQseal Elite Closure Device.

Tray Loading

- A. Whilst the tray is lying flat on the sterile preparation area, fill the recess at the end of the tray (indicated on the tray with "A") with sterile non-heparinised saline solution (approximately 20 cc) to submerge the implant, Figure 2.


Figure 2. Fill tray recess with non-heparinised saline.

- B. Pull off the small Removable tray lid, indicated with a "B" on the tray, by lifting the tab as shown in Figure 3 until completely removed.

Note: do not remove any other lids.

Figure 3. Remove Removable Tray Lid.

- C. Holding the tray steady, grip the rear of the handle (release cam) and pull until a firm stop is reached (indicated on the tray with “C”), Figure 4. The implant is now folded in the loading cannula tube.

Figure 4. Pull device handle back to firm stop.

- D. Flip the funnel lever to a vertical position to release the delivery device from the loading funnel (indicated on the tray with “D”), Figure 5

Figure 5. Flip the funnel lever to release device from tray.

- E. Lifting the Handle End upwards to overcome the stop, continue to pull back and withdraw the PerQseal Elite Closure Device from the tray.
F. Visually check the integrity of the implant (within the Loading-cannula); there should be no creases, folds or tears, see Figure 6.

Figure 6. Inspection of Implant (within loading cannula).

- G. The PerQseal Elite Closure Device is now ready for use.

Note: Use the PerQseal Elite Closure Device within 30 minutes of completing this step.

Procedural Steps:

Inserting the Device:

- a. Use non-heparinised saline to lubricate the guidewire during insertion. Backload the *in situ* guidewire into the device via the implant within the Loading-cannula, Figure 7, continue to insert the guidewire until it emerges from the rear of the PerQseal Elite Closure Device handle. Hold the PerQseal Elite Introducer-sheath at a low angle relative to the longitudinal axis of the artery (e.g. 20°) and keep the Introducer-sheath as straight as practical, to avoid kinking.


Figure 7. Loading guidewire into device cannula.

First Lock:

- b. Advance the Loading-cannula over the guidewire towards the Introducer-sheath hub. Insert the Loading-cannula into the hub of the Introducer-sheath until fully engaged. If difficulty is experienced tracking over the guidewire, ensure the guidewire is not kinked and that it is free of coagulated blood and thrombus. The loading-cannula cap should be fully seated and locked into the Introducer-sheath hub, (with a “click”). The green marking being aligned confirms this step is successfully completed, Figure 8.

Note 1: Hold the guidewire at the rear of the handle during insertion to ensure the device travels over the wire and do not inadvertently insert additional wire into the patient during device insertion/advancement.

Note 2: Only grip the back portion of the Loading cannula during insertion as indicated in Figure 8.



Figure 8. Advance the PerQseal Elite Closure Device to the Introducer-sheath.

- c. Maintaining the position of the Introducer-sheath and guidewire, advance the device handle to engage with the Loading-cannula cap, this will lock the Loading-cannula cap to the handle (with a “click”), the delivery system and introducer sheath are now one unit, Figure 9. The implant is now positioned within the Introducer-sheath at the tip.


Figure 9. Advance PerQseal Elite Closure Device handle to engage with Loading-cannula cap.

Device Actuation:

STEP 1 *Unsheathing the implant*

- a. Rotate the Sheath-cam 360° clockwise (to its non-return stop with a “click”) to withdraw the sheath and expose the implant within the artery, Figure 10.


Figure 10. STEP 1 of PerQseal Elite Closure Device Actuation (Unsheathing the implant).

The green mark on the Loading-cannula cap and Introducer-sheath hub will appear in the window on the front of the handle, Figure 11.

Figure 11. Introducer-sheath fully retracted into handle with green marking within handle window.

STEP 2 *Tamponade*

- a. While holding the device handle at its lowest possible angle (≤ 30 degrees), withdraw the delivery system and sheath in a coaxial orientation relative to the vessel (whilst applying light extracorporeal pressure above the access site) until the implant is adjacent to the arteriotomy in the tamponade position. At this point, gentle tactile resistance should be felt as a result of the implant anchoring itself against the inside lumen of the artery wall. This is the tamponade position, and the implant is now controlling the bleeding from the arteriotomy, Figure 12.


Figure 12. STEP 2 of PerQseal Elite Closure Device Actuation (Tamponade), showing position of implant juxtaposition to the arteriotomy.

Note 1: It is important to confirm the location by gentle resistance of the implant against the arterial wall to avoid the implant being deployed within the tissue tract.

Note 2: Avoid excessive tension on the device. Excessive tension on the device will reduce the implant's effectiveness to achieve haemostasis or could result in the implant being pulled out of the arteriotomy.

- b. Slight adjustment of PerQseal Elite Closure Device position may be necessary to achieve tamponade. Ensure that an acceptable tamponade is achieved before proceeding with deployment.
- c. Open the sidearm three-way tap to confirm tamponade (minimal blood should flow from the sidearm when an acceptable tamponade is achieved).

NOTE: If tension is needed to achieve initial tamponade this should be accompanied with counter-tension at the skin level. If the implant is positioned correctly against the lumen, tamponade should be achieved with minimal tension on the device. The position of the implant relative to the vessel wall/arteriotomy maybe confirmed by ultrasound visualisation.

STEP 3 Depress Cam-lock

Depress the Cam Lock to release the fixation cam, Figure 13.

Figure 13. Step 3 Depress Cam lock

STEP 4 Fixation

- Raise the angle of the device relative to the skin to approximately 45° – 50°, Ensure the angle of the device (relative to the parasagittal plane) follows the artery or sheath in its natural lie.
- With minimal tension on the device, rotate the Fixation Cam 360° clockwise (to its non-return stop with a “click”) to deploy the External Fixation, Figure 14. Avoid rotation of the handle with the distal grip.



Figure 14: STEP 4 of PerQseal Elite Closure Device Actuation (Fixation).

STEP 5 Guidewire Removal

Note: Ensure that haemostasis is acceptable before proceeding with guidewire removal.

Remove the guidewire completely from the patient and PerQseal Elite Closure Device. Guidewire removal may be done by an assistant, Figure 15.


Figure 15. Step 5 Removal of guidewire.

STEP 6 Implant Deployment and Release

- Without delay, once the guidewire is removed, rotate Release Cam 180° clockwise to push home the guidewire pin and release the Implant from the Delivery System, Figure 16. Avoid rotation of the handle with the distal grip.
- Gently withdraw the PerQseal Elite Closure Device Delivery System & the Introducer-sheath from the tissue tract, Figure 17.

Note: Ensure the hand gripping the distal end of the handle remains fixed and the handle remains stationary without rotation during STEPS 4, 5 & 6.

Figure 16. Step 6 Rotate Release cam.

Figure 17. Removal of deployed device.

- If suitable access is available, an angiogram or colour Doppler duplex ultrasound may be performed immediately after PerQseal Elite Closure Device deployment to confirm arteriotomy sealing.

Note: If haemostasis has not occurred, gentle manual compression may be applied proximal to the arteriotomy sufficient to control bleeding. Consult troubleshooting section.

- Use standard clinical procedures to treat any skin incisions.

Post procedural Management

Provide patient with Patient Implant Card containing all of the required information.

Before considering patient discharge, assess for the following:

- Bleeding from the access site
- Anticoagulation thrombolytic or antiplatelet therapy
- Haematoma at the access site
- Hypotension
- Any condition requiring clinical observation
- Pain, discomfort or transitory local irritation and inflammation at the puncture site.
- Local and distal pulse deficits, ischaemia.

Trouble Shooting Guide

1. **Cannot insert the loading cannula into the Introducer-sheath hub:**
Ensure that the supplied PerQseal Elite Introducer-sheath is being used and the loading cannula is not damaged.
2. **Difficult to advance the PerQseal Elite Closure Device into the Introducer-sheath:**
Ensure that the loading cannula is fully engaged with the Introducer-sheath hub, if not discard the PerQseal Elite Closure Device and replace with new PerQseal Elite Closure Device.
3. **PerQseal Elite Closure Device implant withdrawn from artery during positioning:**
Control bleeding by application of manual compression or inflate an endovascular balloon if *in situ*. Withdraw the PerQseal Elite Closure Device and Introducer-sheath from the patient, as per Section 1 below (Recommended Recovery Options), ensuring to leave the guidewire *in situ*. Discard the PerQseal Elite Closure Device and Introducer-sheath and replace with new PerQseal Elite Closure Device.
4. **Cannot achieve acceptable tamponade during implant positioning prior to deployment:**
 - Avoid excessive tension on the PerQseal Elite Closure Device as excessive tension will reduce the implant's effectiveness to achieve haemostasis and control bleeding.
 - Change the angle and/or orientation of the PerQseal Elite Closure Device (relative to the longitudinal axis of the artery) by using very small adjustments (approximately 5°).
 - If acceptable tamponade cannot be achieved do not deploy the PerQseal Elite Closure Device. Follow the recovery suggestion outlined in Section 1 or 2 below (Recommended Recovery Options).
5. **Difficult to withdraw the guidewire from the PerQseal Elite Closure Device:**
 - Do not force.
 - Ensure the guidewire is not snagged endovascularly.
 - Lower the angle of the PerQseal Elite Closure Device (relative to the longitudinal axis of the artery).
 - If it is not possible to withdraw the guidewire do not deploy the PerQseal Elite Closure Device and follow the steps outlined in Section 1 or 2 Recommended Recovery Options below.
6. **Haemostasis not achieved post deployment of the PerQseal Elite Closure Device:**
 - Light manual compression may be applied proximal to the implant to control bleeding until haemostasis is achieved.
 - Note: It is not recommended to apply manual pressure directly onto the location where implant is situated. Excessive compression may damage the implant and adversely cause bleeding.
 - Ensure ACT level is below 350 seconds, use clinical judgement and reduce the effects of heparin with agents such as protamine sulphate.
 - Coagulating extraarterial agents may be used.
 - Temporary use of an endovascular balloon inflated directly beneath the implant may be used.
 - Keep patient under observation until stable.
7. **Difficult to withdraw the delivery device post deployment:**
The device tip may be caught on subcuticular soft tissue. Rotate the device in either direction and try to withdraw.
8. **Guidewire in situ post deployment:**
If the guidewire is still *in situ* post deployment, follow the recovery steps outlined in Section 5 Recommended Recovery Options.

Recommended Recovery Options

1. If it becomes necessary to remove the implant from the artery whilst still attached to the delivery system and prior to the external fixation deployment (Step 4 above) then follow these recommended recovery steps:
 - i) Control bleeding by manual compression or via endovascular balloon inflation proximal to the access site.
 - ii) Ensure Introducer-sheath is withdrawn from the artery (bring device to Step 2 Achieve Tamponade stage).
 - iii) Advance the PerQseal Elite Closure Device into the artery approximately 10 mm.
 - iv) Rotate the handle 180° **anticlockwise** so that the handle is now facing towards the patient in a downward direction.
 - v) Reduce the device angle to < 20° (or as low as possible) relative to the patient, withdraw the handle to remove the PerQseal Elite Closure Device and implant from the artery. Ensure the 0.035" guidewire remains *in situ* during withdrawal.
 - vi) Inspect the PerQseal Elite implant to ensure the scaffold and patch have been removed from the patient.
 - vii) Open a new PerQseal Elite Closure Device procedure pack.
 - viii) Insert a new Introducer-sheath as per PerQseal Elite Introducer IFU.
 - ix) Follow the PerQseal Elite Closure Device deployment steps outlined above.
2. If it becomes necessary to remove the implant from the artery whilst still attached to the delivery system post the external fixation deployment (Step 4 above) then follow these recommended recovery steps:
 - x) Control bleeding by manual compression or via endovascular balloon inflation proximal to the access site.
 - xi) Ensure Introducer-sheath is withdrawn from the artery (bring device to Step 2 Achieve Tamponade stage).
 - xii) Advance the PerQseal Elite Closure Device towards the artery and in a cranial direction relative to the arteriotomy as far as possible, approximately 10 – 15 mm.
 - xiii) Whilst maintaining pressure on the device in the cranial direction, apply digital pressure proximal to the relative position of the implant to collapse the proximal artery and elevate the device handle and shaft to a vertical position.
 - xiv) Maintain the cranial pressure on the device and reduce the device angle to ~30° relative to the patient, withdraw the handle to remove the PerQseal Elite Closure Device and implant from the artery. Ensure the 0.035" guidewire remains *in situ* during withdrawal.
 - xv) Inspect the PerQseal Elite implant to ensure the scaffold and patch have been removed from the patient.
 - xvi) Open a new PerQseal Elite Closure Device procedure pack.
 - xvii) Insert a new Introducer-sheath as per PerQseal Elite Introducer IFU.
 - xviii) Follow the PerQseal Elite Closure Device deployment steps outlined above.
3. If the implant becomes detached from the delivery system within the artery whilst the Introducer-sheath and guidewire are *in situ* use the following recovery steps to retrieve the implant:
 - i) Withdraw the delivery system and the Introducer-sheath and ensure the guidewire remains *in situ*.
 - ii) Re-cannulate the artery over the 0.035" guidewire
 - iii) Open a new PerQseal Elite Closure Device Procedure Pack and insert Introducer as per PerQseal Elite Introducer IFU.
 - iv) Remove the dilator.
 - v) Insert an additional 0.035" guidewire to retain access.
 - vi) Snare the original 0.035" guidewire beyond the implant.
 - vii) Withdraw the snare, original 0.035" guidewire and trapped implant into the Introducer-sheath.
 - viii) Ensure the implant is within the Introducer-sheath and remove the snare, original guidewire, trapped implant and Introducer-sheath from the patient, whilst maintaining the position of the second 0.035" guidewire *in situ*.
 - ix) Insert PerQseal Elite Introducer and deploy the new PerQseal Elite Closure Device as per instructions above.

4. For recovery of the bioabsorbable implant from within the artery post Introducer-sheath and guidewire removal it is recommended to locate the implant via angiography. Based on the implant's location and the level of perfusion downstream of the implant, use clinical judgement to decide an appropriate strategy for recovery.
5. If the device has been fully deployed whilst the 0.035" guidewire is still *in situ* use the following recovery steps to retrieve the guidewire.
 - i) Remove the sheath and delivery system from the guidewire.
 - ii) Very slowly increase tension on the guidewire, at the same angle as deployment or approx. 30° to the artery, until movement of the guidewire occurs. Simultaneously apply light extracorporeal pressure slightly distal to the position of the implant.
 - iii) Continue to withdraw the guidewire whilst maintaining light extracorporeal pressure until fully removed.
 - iv) Inspect the guidewire to ensure the entire guidewire has been removed from the patient.
 - v) Expect some bleeding from the guidewire port in the implant. This can be managed by applying sufficient manual compression to control bleeding.

Note: excessive compressing will damage the implant and adversely cause bleeding.
 - vi) Keep the patient under observation until stable.

Adverse Events & Risks

The following are anticipated complications associated with use of the PerQseal Elite closure device:

- Bruising, oozing at the puncture site
- Transient extravasation
- Late bleed
- Retroperitoneal bleeding
- Haematoma
- Pain, discomfort or transitory local irritation and inflammation at the puncture site
- A mild local inflammatory reaction to the implant material (related to the ongoing biodegradation process) within the first 90 days. The breakdown products of the implant material are glyoxylate or glycine which is extracted via urine or as carbon dioxide respectively
- Access-site related nerve injury, vascular spasm
- Local and/or distal pulse deficits, ischemia
- Stenosis at the closure site
- Infection
- Arterial or deep vein thrombosis
- Pseudoaneurysm
- Arteriovenous fistula
- Vasovagal response

Other events that could possibly occur include:

- Occlusive intraluminal thrombus and/or emboli formation at the implant site
- Sinus formation
- Compartment syndrome
- Iliac or femoral artery dissection
- Partial or complete occlusion of the artery or peripheral arteries
- Embolisation of the intraarterial implant components
- Damage of the arterial wall (perforation or rupture)
- Allergic or adverse foreign body reaction
- Corrective intervention due to any of the above complications such as surgery and/or transfusion

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the procedure was performed.

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