PerQseal[®]
Elite

Introducer

MD

For assistance, contact: Vivasure Medical customer Service: Für Hilfestellungen wenden Sie sich bitte an: Vivasure Medical Kundendienst: Pour toute assistance, contactez: Service client de Vivasure Medical : Neem voor ondersteuning contact op met: Klantenservice Vivasure Medical: Hvis du har brug for hjælp, bedes du kontakte: Vivasure Medical Customer Service: Si necesita ayuda, póngase en contacto con Atención al cliente de Vivasure Medical: Για βοήθεια, επικοινωνήστε: Τμήμα εξυπηρέτησης πελατών της Vivasure Medical	Telephone Number (9:00 – 17:00 GMT Mon to Fri) Telefonnummer (9:00 – 17:00 GMT Mo bis Fr) Numéro de téléphone (9h00-17h00 GMT, du lundi au vendredi) Telefoonnummer (9:00 – 17:00 GMT, ma t/m vr) Katalognummer, (9:00 – 17:00 GMT mandag til fredag) Número de teléfono (de lunes a viernes de 9:00 a 17:00 GMT) Αριθμός τηλεφώνου (Δευτέρα έως Παρασκευή, 9:00 – 17:00 GMT)	+353 91 395 440
Basic UDI-DI, Basic UDI-DI, IUD de base, Basic UDI (identificatiecode), Grundlæggende UDI, UDI básico, Βασικό UDI-DI	539153318DP4FA1INTV7	

En: INSTRUCTIONS FOR USE	2
De: GEBRAUCHSANWEISUNG.....	6
Fr : MODE D'EMPLOI	11
NL: GEBRUIKSAANWIJZING	16
Da: BRUGSANVISNING	21
Es: INSTRUCCIONES DE USO	26
El: ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	31

Device Description

The PerQseal[®] Elite Introducer is a single-use, disposable, radiopaque catheter that provides a flexible conduit for the insertion of the PerQseal[®] Elite Closure Device into the femoral vasculature. The PerQseal[®] Elite Introducer consists of an Introducer-sheath (with a hydrophilic coating) and a Dilator with a blood signal port, which accommodates up to a 0.035" (0.89 mm) guidewire. The blood signalling is via a sidearm attached to the sheath-hub. The 3-way tap on the sidearm allows control of the blood signalling and serves as a sheath flush port. Refer to Figure 1 which shows the key features of the device. The PerQseal[®] Elite Introducer is available in 18 F.


Figure 1. PerQseal Elite Introducer

- | | |
|---|--|
| 1. Dilator (top view) | 6. Introducer-sheath engagement recess |
| 2. Dilator tapered tip Introducer-sheath (top view) | 7. Introducer sidearm |
| 3. Dilator engagement tab | 8. Sidearm 3-way tap (Blood signal and flush port) |
| 4. Introducer-sheath (top view) | 9. Dilator blood inlet hole |
| 5. Introducer-sheath hub (with haemostatic valve) | |

Table 1. Product Size

Product Code	Introducer-sheath Dimensions		
	Effective Length	Nominal Outer Diameter	Nominal Inner Diameter
18 F	mm	mm	mm
	150	7.3	6.3
Product Code	Dilator Dimensions		
	Effective Length	Nominal Outer Diameter	Nominal Inner Diameter
18 F	mm	mm	mm
	273	6.5	0.9

Additional Recommended Equipment

- 20 ml syringe
- 0.035" (0.89 mm) guidewire

How Supplied

The PerQseal Elite Introducer is non-pyrogenic and supplied STERILE. The PerQseal Elite Introducer is provided sterile, using an ethylene oxide process. The PerQseal Elite Introducer is supplied in a Tyvek[®] pouch and mounted on a HDPE card. The Tyvek pouch provides the sterile barrier system. All packaging not contained within the Tyvek pouches are non-sterile. Package content includes the Introducer-sheath, Dilator and Instructions for Use. Store in a cool, dry, dark place.

Please refer to the end of this booklet for graphical symbols for medical device labelling.

Vivasure Medical Limited
Parkmore Business Park West
Galway
Ireland, H91 V3KP
www.vivasuremedical.com

¹ Tyvek[®] is a trademark of DuPont.

Intended Purpose

The PerQseal Elite Introducer is intended to provide access to the femoral vasculature, to facilitate the introduction of the PerQseal Elite Closure Device and to minimize blood loss associated with such procedures.

Indications for Use

The PerQseal Introducer is intended to provide access to the femoral vasculature and facilitate the introduction of the PerQseal Elite Closure Device.

Intended Patient Population

- Individuals over 18 years of age.
- Females who are not pregnant or lactating.
- Patients that are not vulnerable.

Intended User

The intended users of the PerQseal Elite include, but are not limited to, vascular surgeons, interventional radiologists, interventional cardiologists, cardiothoracic surgeons and angiologists. All users are expected to be medically trained and qualified practitioners.

Clinical Benefits

- Less invasive percutaneous access compared to surgical access
- Minimisation of pain and discomfort compared to surgical access
- Minimisation of secondary interventions to control bleeding
- Percutaneous access leads to minimal procedural steps required to achieve access to vessels compared to surgical access
- Percutaneous access leads to shorter overall procedure time compared to surgical access
- Percutaneous procedure has lower major complications rates than alternative therapy of cut-down access
- Safe and Effective vessel access for subjects treated with anticoagulation therapy, antiplatelet agents, intravenous glycoprotein IIb/IIIa inhibitors, or thrombolytic agents
- Facilitates guidewire access throughout procedures
- Reduction in scarring compared to surgical cut-down

Clinical Results

The Frontier V clinical study was performed to confirm safety and performance of an equivalent device, the PerQseal+ Closure Device and PerQseal L Introducers, to percutaneously close femoral artery punctures and to induce arterial haemostasis in patients undergoing endovascular procedures requiring an arteriotomy up to 26 F. The Frontier V study was a prospective, multicentre, non-randomised study to investigate the safety and performance of the PerQseal+. The study was not blinded prior to, during or post the procedure. Patients were not excluded on the basis of age, race, concomitant therapy or co-existing diseases.

Primary Endpoint (Safety): Incidence of major vascular access site complications related to the PerQseal+ up to 1-month from implantation, is no worse than the performance goal for safety estimated from literature (in an equivalent patient population) for percutaneous vascular closure.

Performance: Assessed by technical success rate for the PerQseal+ at discharge is no worse than success rate estimated from literature for alternative large hole closure devices from an equivalent patient population.

Study Population (7 EU sites)	Subjects N (Closures)		Male, N (%)	Age, years
	91 (103)		63 (69%)	78.2 (±8.1)
Primary Procedure	TAVR		EVAR	
	56		35	
Follow Up	1 month		3 months	
	87		84	
Results	PerQseal+ Related Major Complication Rate, N (%)		Technical Success Rate	Time to haemostasis (minutes) [Mean ± SD, Median]
	1 (1%)		95%	4.0 ± 5.0, 1.0

Conclusion: Results from the Frontier V clinical study have demonstrated the safety and efficacy of the PerQseal+ closure device for the percutaneous closure of large-bore femoral arteriotomies.

Some of the limitations of the Frontier V clinical study include:

- The patient population was selected based on the inclusion/exclusion criteria. These exclusions included a restriction of the Sheath to Femoral Artery Ratio greater than 1.05 and anterior or circumferential calcification at or near the access site.
- There are a relatively small number of subjects with a BMI greater than 30 in the study population
- There was a relatively small number of subjects enrolled in the study with a primary sheath size greater than 19 F

The Elite clinical study was performed to confirm safety and performance of the PerQseal Elite Closure Device and PerQseal Elite Introducer, to percutaneously close femoral vein punctures and to induce haemostasis in patients undergoing endovascular procedures requiring a venotomy up to 26 F. The Elite clinical study was a prospective, multicentre, non-randomised study to investigate the safety and performance of the PerQseal Elite. The study was not blinded prior to, during or post the procedure. Patients were not excluded on the basis of age, race, concomitant therapy or co-existing diseases.

Primary Endpoint (Safety): Incidence of combined of both major and minor vascular access site complications attributable to the PerQseal Elite device through 30 days.

Primary Endpoint (Effectiveness): Time to Haemostasis (TTH) defined as elapsed time in minutes from PerQseal Elite (Introducer sheath and Delivery device) removal to first observed cessation of common femoral vein (CFV) bleeding (clinically acceptable cessation of venous bleeding), excluding cutaneous or subcutaneous oozing, and in the absence of a developing hematoma.

Analysis of primary endpoints

For both primary endpoints, a null and alternative hypothesis were tested. Both null hypotheses need to be rejected for the study to be considered successful (i.e., noninferiority can be claimed). The hypotheses:

- H0 for safety: Rate of device-related complications = 0.187*
- H1 for safety: Rate of device-related complications < 0.187*
- H0 for effectiveness: TTH = 16.97 minutes*
- H1 for effectiveness: TTH < 16.97 minutes*

*Rate estimated from literature for alternative large hole closure devices from an equivalent patient population.

Success on the primary safety and effectiveness endpoint has been established. There were no complications or adverse events associated with the use of the PerQseal Elite Introducer.

Study Population	Subjects N (Closures)	Male, N (%)	Age, years
Centres (Germany & Netherlands)	46 (47)	30 (65.2%)	77.4 ± 8.1
Follow Up	1 month		
	44		
Results	PerQseal Elite Related Complication Rate, N (%)	Technical Success Rate, N (%)	Time to haemostasis (minutes) [Mean(SD), Median (min, max)]
	1 (2.3%)	47/47 (100%)	0.15 (0.55), 0 (0, 3.00)

Conclusion: Results from the PerQseal Elite clinical study have demonstrated the safety and efficacy of the PerQseal Elite closure device for the percutaneous closure of large-bore femoral Venotomies. The study does contain a few limitations. The study was a single-arm study, thereby precluding any direct comparison with other marketed VCDs. The patient population was selected based on the inclusion/exclusion criteria agreed upon by the regulatory body approving the study. As a result, the selected patient population may not fully represent the diversity and clinical complexities of the real-world patient population who require vascular closure. Post Market Clinical Follow-up studies will continue in order to confirm the safety and efficacy of the PerQseal Elite in a larger patient population.

Contraindications

There are no contraindications to the use of this device. Attention is drawn to the warnings and precautions.

Warnings

- Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, call your VIVASURE representative.
- For single use only. Do not reuse, reprocess or resterilise. Reuse, reprocessing or resterilisation may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
- Do not use the Introducer-sheath and/or Dilator if the packaging has been previously opened or damaged, or if any of the components appear to be damaged or defective. This may result in patient injury, illness or death.
- Discard product if Use By Date has passed or Use By Date is illegible.

Precautions

- This product is intended for use by physicians trained and experienced in diagnostic and therapeutic endovascular techniques. Standard techniques for placement of vascular access Introducers should be employed.
- Do not advance, torque or withdraw the Introducer-sheath and Dilator against excessive resistance. Severe vessel damage and/or injury and/or breakage of the device may occur.
- Do not advance the Introducer-sheath without the Dilator fully engaged with the Introducer-sheath. Do not advance the Introducer without an *in situ* 0.035" guidewire. Severe vascular damage and/or injury may occur.
- All devices used with the Introducer-sheath should move freely through the valve and sheath. Damage to the valve/sheath may result if the fit is tight.
- Do not use the PerQseal Elite Introducer where bacterial contamination of the surrounding tissues may have occurred.
- Do not use the PerQseal Elite Introducer if the puncture is through the posterior wall or if there are multiple punctures, since such punctures may result in a retroperitoneal bleeding.
- Do not use the PerQseal Elite Introducer if venous access is or is suspected of being, via an artery.
- Do not use the PerQseal Elite Introducer if there is a vascular graft or stent within 20 mm of the access site.
- Do not use the PerQseal Elite Introducer if the Introducer-sheath to vessel diameter ratio is greater than or equal to 1.05. Do not use PerQseal Elite Introducer if vessel lumen diameter is less than 7 mm.
- Do not use in patients < 18 years of age.
- As with all percutaneous based procedures, infection is a possibility. Observe sterile techniques at all times when using the PerQseal Elite Introducer. Employ appropriate groin management post procedure and post hospital discharge to prevent infection.
- If a patient has had a procedural sheath left in place for longer than 8 hours, consider the use of prophylactic antibiotics before use of the PerQseal® Elite Introducer.
- Consideration should be given to the access site: tortuosity, disease state (including calcification, plaque and thrombus) to ensure successful Introducer-sheath introduction and subsequent withdrawal. If the vessel is not adequate for access, major bleeding, vessel damage, or serious injury to the patient, including death, may result.

PROCEDURE

The PerQseal Elite Introducer is to be used only by a licensed physician possessing documented instruction and training by Vivasure Medical in the use of the PerQseal Elite Introducer.

Examination of Packaging & Labelling

- Inspect visually the sterile barrier system (Tyvek Pouch) to ensure there is no damage or breaches of packaging integrity.
- Check all labelling to ensure the product is within its use by date.

Opening the packaging

- While observing aseptic techniques, open the Tyvek pouch by peeling apart the end with the chevron seal (the end with the 'v-shape' seal).
- Present the exposed card (with Introducer-sheath and Dilator attached) for aseptic transfer into the sterile field.

Preparation

- Flush the Dilator by injecting saline through the Dilator-hub and ensure the tip of the Dilator tapered end (2, Fig. 1) is wet prior to insertion of the Dilator into the Introducer-sheath.
- Insert the Dilator into the Introducer-sheath by inserting the tapered end (2, Fig. 1) in through the Introducer-sheath hub (5, Fig. 1), hold the Dilator shaft as close as practical to the sheath hub during insertion. Ensure the engagement tab (3, Fig. 1) on the Dilator-hub aligns with the engagement recess (6, Fig. 1) in the Introducer-sheath hub. Fully insert the Dilator into the Introducer-sheath to engage the snaps on the Dilator-hub with the sheath hub. The blood inlet hole (9, Fig. 1) on the Dilator should be exposed forward of the Introducer-sheath (4, Fig. 1).

- c. Flush Introducer with saline via the sidearm 3-way tap. Ensure the flush exits the Dilator tip and blood inlet hole (9, Fig. 1).
- d. Activate the hydrophilic coating on the Introducer-sheath by wetting the outer surface of the Introducer-sheath with saline (either heparinised or non-heparinised). Ensure the Introducer-sheath remains wet during placement.
- e. Leave the sidearm 3-way tap in the open position for blood signalling during insertion.
- f. Follow standard clinical practice for vessel puncture or sheath exchange as required.

Insertion into vessel

- a. Ensure the graduated markings on the PerQseal Elite Introducer are facing skyward as the graduated markings are required to attain correct orientation during subsequent procedure steps.
- b. Insert the PerQseal Elite Introducer over the in-situ 0.035" guidewire until blood signal is observed from the side-arm 3-way tap (9, Fig. 1), stop advancing and slowly withdraw the Introducer until the blood signal first stops, stop withdrawal at this point and close the side-arm 3-way tap.
- c. Take a note of the number on the graduated scale of the PerQseal Elite Introducer at the skin level (or vessel wall if vessel exposed) when this blood signal first occurs. This gives an indication when the tip of the Introducer-sheath is in the vessel.
Note: The introducer may be withdrawn and re-advanced to confirm correct position.
- d. At a low angle ($\leq 30^\circ$), advance the PerQseal Elite Introducer a further **4 graduations** into the vessel from the blood signal point. Stop advancing at this point.
- e. Squeeze the dilator side tabs and remove the Dilator while maintaining the 0.035" guidewire and the introducer-sheath in its position within the vessel. Maintain the Introducer in this position until it is locked with the PerQseal Elite Closure Device.
- f. Flush the Introducer-sheath sidearm with non-heparinised saline.
Note: Maintain the Introducer-sheath in-line with the vessel and as straight as possible to prevent kinking of the Introducer-sheath
- g. Follow standard clinical practice for guidewire insertion and exchange (as required).

Adverse Events & Risks

The following are anticipated complications associated with use of the Introducer:

- Allergic response to materials
- Blood loss, bleeding or haematoma
- Retroperitoneal bleeding
- Embolisation (micro or macro) with transient or permanent ischemia
- Infection
- Vascular trauma (e.g. dissection, rupture, perforation, or tear)
- Procedural discomfort
- Thrombosis
- Transitory vessel occlusion

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the procedure was performed.

Warranty

VIVASURE warrants that reasonable care has been used in the design and manufacture of this device. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilisation of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond VIVASURE's control directly affect the device and the results obtained from its use. VIVASURE shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. VIVASURE neither assumes, nor authorises any other person to assume for it, any other or additional liability or responsibility in connection with this device. **VIVASURE assumes no liability with respect to devices reused, reprocessed or resterilised and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such devices.**

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